

Liquid Handling in Perfection

COMFORT AND EFFICIENCY IN SERIAL PIPETTING



VITLAB
Competence in Labware

Pipet controller for the most demanding users.

Pipet controller VITLAB pipeo®

VITLAB specialises in accurate dispensing, titration and pipetting, and offers a variety of precision instruments for Liquid Handling applications.

With the new pipet controller VITLAB pipeo®, pipet handling is simple and comfortable. The ergonomic handle, very light weight – only about 190 grams – and excellent balance all contribute to operating ease. Freely and precisely adjust the speed of the liquid discharge with one hand using the two easily accessible buttons.

A 50 ml pipet can be filled, quietly and with minimum vibration, in less than ten seconds. The VITLAB pipeo® can be set to release liquid either in gravity delivery, or in blow out mode using the battery-operated motor.

To avoid surprises, the LED lamp changes from green to red approx. two hours before the battery needs to be recharged. If necessary, the VITLAB pipeo® can be used as the battery is recharged so that important work is not interrupted.



Charge-level indicator (LED)

Battery compartment cover

Mode selector (gravity delivery / blow-out)

Single-hand operation

Direct discharge of liquid vapors

Adapter for a secure attachment of the pipets, from 0.1 to 100 ml

Pipetting buttons

Item no.	
Pipet controller VITLAB pipeo®	
1631500	VITLAB pipeo® with Battery Charger for continental Europe 230 V/50 Hz
1631510	VITLAB pipeo® with Battery Charger for UK and Ireland 230 V/50 Hz
1631520	VITLAB pipeo® with Battery Charger for Australia 240 V/50 Hz
1631530	VITLAB pipeo® with Battery Charger for Japan 100 V/50 Hz

Spare parts are available upon request.