



Partner per l'Italia



STEDA SRL
Marketing Trading Logistics



Table of Contents

Everything you need for a well stocked laboratory
or educational environment

| PG # | Description | Summario |
|------|-------------------------------|--|
| 4 | MAG FUGE | MagFuge®, Centrifuga & Agitatore |
| 5 | GUSTO HIGH SPEED CENTRIFUGE | Gusto®, Centrifuga, mini |
| 6 | SPROUT PLUS MINI CENTRIFUGE | Sprout®Plus, Centrifuga, mini |
| 7 | PORTABLE CENTRIFUGE KIT | |
| 8 | MAGNETIC INDUCTION STIRRER | Agitatore, magnetic induction |
| 9 | MINI MAGNETIC STIRRER | Agitatore, magnetic mini |
| 10 | MINI VORTEXER | Agitatore, Vortexer mini, scuotimento |
| 11 | VORTEXER | Agitatore, Vortexer per provette |
| 12 | PIPETTE CONTROLLER | Aspiratore, Pipettatore Controller |
| 12 | PIPETTE CONTROLLER RF1000 | Aspiratore, Pipettatore RF1000® |
| 13 | PIPETTE CONTROLLER RF3000 | Aspiratore, Pipettatore RF3000®, semiautomatico |
| 14 | E COUNT, COLONY COUNTER | eCount, Penna contacolonie digitale |
| 15 | MINI LIGHT, MINI-MAGNIFIER | Mini Light box e Mini lente |
| 15 | HAND TALLY COUNTER | |
| 15 | BENCHTOP COUNTER | Contacellule manuale |
| 16 | LABORATORY NOTEBOOK | Quaderni da laboratorio, a righe e quadretti |
| 16 | TWIN TIP PERMANENT PENS | Penna contacolonie> Penna, colore nero |
| 17 | LAB ALERT CLOCKS & TIMERS | Contaminuti display digitali |
| 18 | LAB DRYING RACK | Scolavaterria da parete-banco, con pioli |
| 19 | WORKSTATION LAB MAT | Tappetini in gomma, protettivi superfici |
| 19 | SILICONE LAB MAT | |
| 20 | BLOOD COLLECTION TRAY | Cestello robusto porta oggetti, con manico |
| 20 | TOTEMAX BLOOD COLLECTION TRAY | ToteMax® Cestello porta oggetti, con manico |
| 21 | DROPLET BLOOD COLLECTION TRAY | Cestello porta oggetti, con divisori, per provette |
| 21 | DROPLET LID | |
| 22 | DURAPORTER | DuraPorter®, contenitore trasporto |
| 23 | DURAPORTER XL | DuraPorter® XL, contenitore trasporto |
| 24 | ALUMINIUM WEIGHING BOATS | Capsule in alluminio, per pesate |
| 24 | LEVGO SMARTPANS | |
| 24 | WEIGHING PAPERS | |
| 25 | SMARTBOATS | smartBoats™, Capsule carta ECO plasticata |
| 25 | LEVGO SMARTSPATULA | smartSpatula™, Spatole monouso in PP ECO |
| 26 | WEIGHING BOATS | BoatsLab, Navicelle in PS antistatico, per pesate |
| 28 | 32 PLACE PCR RACK | Portamicroprovette, 32 PCR rack, strips 8x4 |
| 29 | 50/81/100 WELL MICROTUBE | Portamicroprovette criogenia, 50-81-100 posti |
| 30 | REVERSIBLE TUBE RACK | Portamicroprovette, Doppio foro |
| 30 | ROTA-RACK DUO TUBE RACK | Portamicroprovette, Rotante |
| 30 | TWO TIER MICROTUBE RACK | Portamicroprovette, Incastro |
| 31 | MICROTUBE RACK 80-WELL | Portamicroprovette, Impilabile 80-posti |
| 31 | MICROTUBE STORAGE RACKS | |
| 31 | MICROTUBE RACKS | Portamicroprovette, Robusto con manici |

| PG # | Description | Summario |
|------|--|--|
| 31 | STORAGE RACK FOR 80-WELL | BOX acrilico, scaffale portaprovette 8-posti |
| 32 | ROTA RACK | Portaprovette, rotante e cubico |
| 32 | CUBE TUBE RACK | |
| 33 | 4-WAY TUBE RACK | Portaprovette incastro, 4 lati, PP Ø diversi |
| 34 | TRANSFORMER TUBE RACK | Portaprovette, Cubico |
| 34 | ONERACK MULTI | Portaprovette, OneRack® Multi |
| 35 | ADAPT A RACK | Portaprovette, Adapt-A-Rack® |
| 35 | UNIVERSAL TUBE RACK | Portaprovette, Universal |
| 36 | POP-UP TUBE RACK | Portaprovette, Pop-Up |
| 36 | FOLD & SNAP TUBE RACK | Portaprovette, Fold & Snap |
| 37 | ONERACK | Portaprovette, OneRack® |
| 38 | MEGA RACK | Portaprovette, Mega Rack® |
| 38 | PEG TUBE RACK | Portaprovette, Pioli |
| 39 | TRACRACK | Portaprovette, TracRack® tube Rack |
| 39 | LOW-TEMP PCR RACK | Portaprovette, Low temperature con coperchio |
| 40 | 50/81/100 WELL MICROTUBE STORAGE BOXES | Scatole portamicroprovette, criogenia 50-81-100 posti |
| 41 | CRYOGENIC VIAL BOXES | Scatole criogeniche, in cartone idrorepellente |
| 41 | MICROTUBE STORAGE BOXES | |
| 42 | 100 PLACE PREMIUM PLUS SLIDE BOX | Scatole portavetrini, Premium Plus, 100 posti |
| 43 | PREMIUM MICROSCOPE SLIDE BOXES | Scatole portavetrini ABS, 25-50-100 posti |
| 44 | MICROSCOPE SLIDE BOX STORAGE RACK | |
| 44 | STANDARD MICROSCOPE SLIDE BOXES | Scatole portavetrini Standard e Low temperature |
| 44 | LOW TEMP MICROSCOPE SLIDE BOXES | Box acrilico, Scaffale per scatole portavetrini |
| 45 | CARDBOARD SLIDE TRAY | Portavetrini a libro in cartone pressato |
| 45 | SLIDE TRAY CLASPS | |
| 45 | SLIDE MAILERS | Portavetrini cartone, posti 1-2-4 |
| 46 | COOL CONTAINERS | Criocontenitori True North®, con coperchio |
| 47 | LEVGO ECO SMARTFUNNEL | Imbuti monouso, carta biodegradabile |
| 47 | LEAD RINGS | Anello blocca beute-matracci |
| 47 | SAFETY BOTTLE CARRIERS | Contenitori di sicurezza, porta bottiglie |
| 48 | BIOHAZARD DISPOSAL BAGS | Sacchetti monouso PP, per sterilizzazione |
| 48 | AUTOCLAVABLE DISPOSAL BAGS | |
| 48 | SHARP CONTAINERS | Contenitori per lamette, siringhe, aghi |
| 49 | PPE DISPENSER & HOLDERS | Dispenser porta cuffie, tappini, guanti, occhiali |
| 50 | PPE DISPENSER & HOLDERS | Contenitore acrilico porta dispenser fazzoletti, puntali |
| 51 | GLOVE BOX HOLDERS | Dispensatore acrilico trasparente per scatole guanti |
| 52 | HEAT SEALER WITH CUTTER | Termoasaldatrice elettrica con taglierina |
| 52 | PARAFILM | Parafilm, pellicola sigillante |
| 53 | PARAFILM DISPENSERS | Dispenser per nastro Parafilm |

MagFuge®

A high-speed mini centrifuge and magnetic stirrer in one unit



Holds
5.0 mL
Tubes



EASY SWITCH OVER
NO
TOOLS REQUIRED

High-Speed Mini Centrifuge

Two Units
in One!

Magnetic Stirrer

High-speed Centrifuge is easy to use and designed to deliver maximum speed (12,500rpm/9,783x g) and flexibility with rapid acceleration and fast controlled breaking

Ideal for most protocols requiring fast spins including PCR, micro- filter cell separation and HPLC. Versatility to handle 0.2, 0.5, 1.5, 2.0 and 5.0 mL tube sizes

Intuitive One-Dial Programming

- Set speed and time with ease:
- Adjustable speed from 500 to 12,500 rpm
- Speed setting for RPM or RCF
- Adjustable time from 30 seconds to 60 minutes

Safe to Use

- Safety features help ensure it will not operate when the lid is open, a rotor is jammed, or the unit is out of balance

High-performance Magnetic Stirrer is designed using high quality rare earth magnets creating a robust magnetic coupling that significantly reduce the chance of spin-outs

Perfect for work with sensitive samples and protocols up to 3L requiring aggressive mixing for full dilution

Easy to Use

- Set stir mode, speed and time with easy one-dial programming:
- Select Mode: clockwise, c. clockwise and oscillate
- Adjustable speed from 50 to 2,500 rpm
- Continuous and adjustable time settings from 30 seconds to 60 minutes

Durable Accessories Included

- Re-usable, easy to clean silicone mats included to keep vessel stable and prevent spills and two NdFeB stir bars (40 x 12 mm) with chemical resistant PTFE coating included

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|-------------|-----------------|------------------|
| LN-72045244 | Grey/Purple | 9.2 x 7.3 x 4.7 | 23.5 x 18.5 x 12 |
| LN-72045245 | Grey/Blue | | |

Replacement parts

| | |
|-------------|---|
| LN-72045246 | MagFuge® Rotor, 12 x 1.5/2.0 mL tube |
| LN-72045247 | MagFuge® Rotor, 6 x 5.0 mL tube |
| LN-72045253 | MagFuge® Stirrer Rotor |
| LN-72045254 | MagFuge® Silicone Mat |
| LN-72045260 | Power Adapter |
| LN-72045150 | Tube Adapters, 0.2 mL, 0.5 mL, Pack of 12 (6 ea.) |
| LN-72045148 | Universal Power Cord Set |

Patent Information

Community Design No. 003515337-0001-0002


Euro. Pat. App. 3 246 088

US Design Patent No. D814,045

US Patent No. 10471439

**Package Includes**

| | | |
|--|---|----------------------------------|
| Centrifuge/Stirrer | Stirrer Rotor | 12 Tube Adapters 0.2 mL & 0.5 mL |
| 5.0 mL Tube Rotor | 2 Stir Bars, 35 mm x 12 mm, cylindrical shape | |
| 1.5/2.0 mL Tube Rotor | 2 Silicone Mats for Stirrer Lid | |
| 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs | | |

**Specifications: Centrifuge**

| | |
|------------------------|---|
| Speed Range | 500 – 999 rpm (increments of 1); 1,000 – 5,000 rpm (increments of 5); 5,000 – 12,500 rpm (increments of 25) |
| RCF | 16 – 9,783 x g (max RCF calculated) |
| Run Time Ranges | 30 seconds – 60 minutes (5 second increments) |

Specifications: Magnetic Stirrer

| | |
|-----------------------|---|
| Speed Range | 50 – 999 rpm (increments of 1); 1,000 – 2,500 rpm (increments of 5) |
| Run Time Range | Stirrer Mode 30 seconds to 60 minutes or continuous (5 second increments) |
| | Oscillating Stirrer Mode 1 second – 60 minutes or continuous |
| | Oscillating Period 30 seconds – 15 minutes |
| Capacity | 3 L plus, Diameter of vessel not to exceed 130 mm |



5
Year
Warranty

Back-lit Display

showing time, speed and operational information on the same screen

EASY SWITCH OVER
NO
TOOLS REQUIRED

Adjustable
time and speed

Quick-spin Button

for momentary spin requirements, press and hold to operate, release to end spin cycle

Gusto® High-Speed Mini Centrifuge

Maximum speed with minimal space

This compact centrifuge is ideal for protocols calling for fast spins with a higher RCF requirement

Made of durable ABS with UV-Resistant coating

User Friendly

- Intuitive Push button programming and bright LCD screen make it easy to use

Quiet Operation

- Brush-less motor enables quiet operation and low vibration

Safety Features

- Lid safety feature prevents opening while rotor is turning and clear polycarbonate allows user to see rotor movement for added visual safety
- Includes automatic shut-off feature for sample imbalance and tilt

Specifications

| | |
|------------------------|-------------------------------------|
| Speed Range | 500–12,500 rpm |
| Run Time Ranges | 1 to 30 minutes |
| RCF | 15 - 9,783 x g (max RCF calculated) |




| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--------|-----------------|--------------------|
| LN-72045602 | Blue ■ | 8 x 6.7 x 4.5 | 20.3 x 17.1 x 11.4 |

Replacement parts

| | |
|-------------|--|
| LN-72045148 | Universal Power Cord Set |
| LN-72045680 | Gusto® Rotor 12 x 1.5/2.0 mL Tube (with Cover and Knob) |
| LN-72045682 | Gusto® Rotor 4 x 8 Strip 0.2 mL Tube (with Knob) |
| LN-72045150 | Tube Adapters, 0.2, 0.5 mL, Pack of 12 (6 ea.) |
| LN-72045689 | Gusto® Rotor Cover (for 12 x 1.5/2.0 mL Tube Rotor Only) |
| LN-72045260 | Power Adapter |

Package Includes

Gusto® High-Speed Mini Centrifuge
 12-Place 1.5/2.0 mL Tube Rotor
 32-Place 0.2 mL Tube Rotor (singles, or 4 x 8 Strip)
 12 Tube Adapters 0.2 mL
 12 Tube Adapters 0.5 mL
 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



The lid can be unlocked and opened with one hand



5
Year
Warranty

EASY SWITCH OVER
NO
TOOLS REQUIRED



Sprout® Plus Mini Centrifuge

Mini but Mighty - ideal for small runs

Easy to Use

- Tool-free SnapSpin rotors and quick-spin options make the centrifuges fast and easy to use

Ergonomic Design

- The lid was functionally designed to fit in the palm of the hand for one handed lock and unlock allowing for quick tube extraction
- Large lid size accommodates longer and skirted tubes and tall tube lids
- The housing is made from ABS plastic and lid is polycarbonate plastic enabling the unit to be easily cleaned

Safety Measures

- Fitted with a braking function, the unit quickly stops spinning when powered off or when lid is opened
- A weighted, rubberized base eliminates potential movement and dampens vibration and noise



Specifications

| | |
|--------------------|----------------------------------|
| Speed Range | Approx. 6,000 rpm (fully loaded) |
| RCF | 2,000 x g |



Package Includes

| |
|--|
| Sprout® Plus Mini Centrifuge |
| 6-Place 1.5/2.0 mL Tube Rotor |
| 16-Place 0.2 mL Tube Rotor (singles, or 2 x 8 strip) |
| 6 Tube Adapters 0.2 mL |
| 6 Tube Adapters 0.5 mL |
| 4 Lid O-Rings |
| Storage Case for Rotors and Adapters |
| Tube Rack |
| 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs |



Sprout® Plus

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--|-----------------|--------------------|
| LN-72045279 | Green ■ | 6 x 5 x 4.3 | 15.3 x 12.8 x 10.9 |
| LN-72045280 | Blue ■ | | |

Replacement parts

| | |
|-------------|---|
| LN-72045692 | Sprout® Plus Rotor SnapSpin 6 x 1.5/2.0 mL Tube |
| LN-72045694 | Sprout® Plus Rotor SnapSpin 2 x 8 Strip, 0.2 mL |
| LN-72045693 | Sprout® Plus Lid O-Ring, Pack of 10 |
| LN-72045188 | Sprout® Plus Power Adapter |
| LN-72045150 | Tube Adapters, 0.2 mL, 0.5 mL, Pack of 12 (6 ea.) |



Specially Designed Inserts

seat the centrifuge, all its accessories, and includes a tube rack

12-Volt Car Adapter
for on-the-go testing



Portable Centrifuge Kit

Spin samples *on-the-go*

The Portable Centrifuge Kit offers a convenient way to spin down and prepare samples outside of a conventional laboratory environment

Immediate Field Processing

- Process samples in remote areas before deterioration can take place
- Allows for immediate on-site testing for faster results
- Made of durable ABS plastic for easy clean in field environments

Easy to Use

- Centrifuge is equipped with tool-less SnapSpin rotors for easy change out
- Simple operation with lock/unlock activation

Portable Convenience

- All items in kit are housed in a handled DuraPorter® Transport Box for portability
- The handled box has an air-tight seal for optimal sample protection

See page 22 to learn more about the
DuraPorter® Transport Box

Item No.

LN-72045179 Portable Centrifugation Kit, 100-240 VAC, 50/60Hz Universal Plug + 12V Car Adapter

*Adapter can be used with Sprout® Plus Mini Centrifuge



Package Includes

Sprout® Plus Mini-Centrifuge

6-Place 1.5/2.0 mL Tube Rotor

16-Place 0.2 mL Tube Rotor (singles, or 2 x 8 strip)

6 Tube Adapters 0.2 mL

6 Tube Adapters 0.5 mL

4 Lid O-Rings

Storage Case for Rotors and Adapters

Tube Rack

1 Low Voltage, Double Insulated Power Adapter
(100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



DuraPorter® Transport Box

Tube Rack Insert
(12-Place 1.5/2.0 mL, 14-Place 0.2 mL)

Transport Box Dividers

Sprout® Plus Cradle

Portable 12V Sprout® Plus Centrifuge Car Adapter

Magnetic Induction Stirrer

Strong, motor-less and maintenance free

Superior Functional Design

- Advanced magnetic coil technology drives the spinning of the stir bar maintaining continuous, constant speed for reproducible results
- Magnetic coupling is optimized by using alternating current (AC) as its driving force strong enough to stir viscous solutions
- Motor-less operation generates minimal heat ideal for stirring temperature sensitive samples and environments for long periods
- With no moving parts, produces less vibration and quieter operation over traditional magnetic stirrers

Intuitive Operation

- Multiple stirring modes are easy to program
- Bright digital display precisely controls time and speed
- Quickly meets stirring speeds with rapid acceleration and braking (less than 10 seconds)



Protective silicone anti-skid plate cover for improved safety



Durable ABS housing is corrosion resistant and easy to clean



Specifications

| | |
|-----------------------|---|
| Speed Range | 50–2,000 rpm variable, 1 rpm increments |
| Run Time Range | 30 seconds to 60 minutes, variable, 1 second increments or continuous |
| Stirrer Modes | Single clockwise Single counter clockwise Auto reverse rotation 1 to 99 cycles 50% duty |
| Capacity | 3 L plus |
| Mixing Surface | 165 mm diameter, silicone mat 160 mm diameter |



Package Includes

- Rare Earth, NdFeB Magnetic Induction Stirrer
- 2 Stir Bars 35 mm x 12 mm Cylindrical Shape
- 1 Silicone Anti-Skid Plate Cover
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|-------------|-----------------|-------------------|
| LN-72045255 | Grey/Purple | 9.7 x 7.7 x 3 | 24.5 x 19.5 x 7.5 |
| LN-72045268 | Grey/Blue | | |

Replacement parts

| | |
|-------------|------------------------------------|
| LN-72045269 | MI Stirrer Silicone Mat |
| LN-72045262 | MI Stirrer Universal Power Adapter |



Mini Magnetic Stirrer

Affordable solution, ultimate function

The Mini Magnetic Stirrer is an affordable, dependable solution ideal for general research, academic, life sciences, industrial and clinical lab settings

Durable Design

- Constructed of chemical resistant ABS and PC blend
- Compact footprint saves valuable space on benchtop
- Lightweight for easy transport

Ultimate Function

- Powerful magnetic coupling ensures accurate alignment, reducing spin-outs during acceleration or deceleration
- Performs routine stirring applications with low vibration and quiet operation

U.S. and E.U. Utility Patent Pending
U.S. Design Patent No. D814,045

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|---|------------------|-------------------|
| LN-72045224 | Grey/Purple  | 7.2 x 5.2 x 1.85 | 18.3 x 13.2 x 4.7 |
| LN-72045225 | Grey/Blue  | | |

Replacement Part

LN-72045133 Mini Stirrer Universal Power Adapter

Specifications

Speed Range 200–2,200 rpm (Infinitely adjustable)

Capacity Up to 3 L




Package Includes

Mini Magnetic Stirrer

2 Stir Bars 35 mm x 9 mm, Octagonal Shape

1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs





Mini Vortexer

Small footprint with maximum power

Simple, touch start operation offers exceptional mixing performance with one hand

Robust Design

- Low profile design and small footprint saves valuable space
- Made of ABS and silicone materials that resist most commonly used chemicals
- Stable platform and elastomer feet ensure unit stays stable during operation

Powerful Function

- Rapid vortex action with low activation force
- Quiet operation and low vibration



U.S. and E.U. Utility Patent Pending
U.S. Design Patent No. D814,045

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|-------------|-----------------|-----------------|
| LN-72045288 | Grey/Purple | 5.1 x 4.3 x 3 | 13 x 11 x 7.5 |
| LN-72045289 | Grey/Blue | | |
| LN-72045194 | Black/Black | | |

Replacement parts

| | |
|-------------|--|
| LN-72045195 | Mini Vortexer Silicone Cover, set of 2 |
| LN-72045196 | Mini Vortexer Universal Power Adapter |


Specifications

| | |
|----------|------------------|
| Speed | 4,500 rpm |
| Orbit | 6 mm |
| Capacity | Up to 50 mL tube |



Package Includes

- Mini Vortexer
- 2 Silicone Caps
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs





Vortexer

Multi-purpose, maximum convenience

Innovative Design

- One mixing head does it all with this sleek, modern and compact designed vortex mixer
- Hold slot configurations for multiple sizes of tubes, vessels and assay plates
- Center pad for quick hand-held mixing

Easy to Use

- Control switch with three operation selections (intermittent, off, continuous)
- Variable Speed Rotary Dial provides pre-set speed control selection during shaking mode

Durable Construction

- ABS material resist most commonly used chemicals and is easy to clean
- Stable elastomer feet dampen vibration and prevent movement during operation



Center pad area for hand held mixing



Assay plate x 1
(to be used with retention cords)



25 mL, 50 mL conical tubes x 2
(horizontally mounted)



15 mL conical tubes x 2
(horizontally mounted)



5 mL tubes x 2
(horizontally mounted in the same retention space as 15 mL)



0.5 mL microfuge tubes x 6
(vertically mounted)



1.5/2.0 mL microfuge tubes x 4
(vertically mounted)



0.2 mL microfuge tubes
(strip tubes) x 24
(vertically mounted)

US Design Patent No. D787,085
US Utility Patent No. 9895670

| Item No. | Speed Control | Power Input | L x W x H in | L x W x H cm |
|-------------|---------------|------------------------------|--|---|
| LN-72045213 | Infinity Dial | 220-240 VAC, 50 Hz, 44 W, CE | 6.7 x 7.5 x 7.8 (with mixing head in place) | 17.1 x 19 x 19.7 (with mixing head in place) |

Replacement parts

| | |
|-------------|--------------------------------|
| LN-72045199 | Vortexer Multi-Head |
| LN-72045200 | Vortexer Central Pad, 3 pk |
| LN-72045201 | Vortexer Retention Cords, 3 pk |

Specifications:

| | |
|-------------|--|
| Speed Range | Up to 3,000 rpm |
| Orbit | 3.7 mm |
| Control | 3 Position Control — (top) intermittent, (middle) off, (bottom) continuous |

Pipette Controller

Precise fluid control at your finger tips

Comfortable to hold and easy to use, the manual inline pipette controller was designed with an ergonomic grip and bulb for simple and precise operation

A T-bar thumb lever controls the aspirating and dispensing speeds from serological pipettes from 1 to 100 mL in size

- Conveniently accommodates either glass or plastic serological or volumetric pipettes ranging from 1–100 mL
- Includes a replaceable .45 µm hydrophobic filter to protect your samples and internal components



Nose cone assembly only

| Item No. | Color | Volume Range mL | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|---------------|----------------|
| LN-72045239 | Blue | 1–100 | 9 x 2.1 x 2.1 | 23 x 5.4 x 5.4 |

Operation Is Simple:

Squeeze the bulb

to create a vacuum to draw up fluid

Push the control lever forward (down)

to dispense the liquid

Push the control lever backward (up)

to aspirate liquid through a pipette

Press the central bulb on the unit

to blow out any remaining fluid if necessary



1
Year
Warranty

Filter not included
in warranty

Stay Organized

with three additional nose cones (green, white, blue) for color identification for both the Pipette Controller & RF1000



Squeeze
to create a vacuum



Filter not included
in warranty

RF1000

Lightweight and easy to use

The RF1000 manual pipette controller is uniquely designed for one-handed operation

Simply squeeze the main bulb and push the toggle forward to aspirate or back to dispense. Push the blowout bulb to expel any liquid remaining in the pipette

- Tapered silicone rubber pipette adapter allows a large range of pipette sizes ranging 1–100 mL and is replaceable
- Includes a replaceable .45 µm hydrophobic filter to protect your samples and internal components



Nose cone assembly only

| Item No. | Color | Volume Range mL | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|---------------|
| LN-72045497 | Blue | 1–100 | 7.1 x 1.8 x 4.7 | 18 x 4.5 x 12 |

Replacement Parts

| Item No. | Description | Color |
|-------------|--------------------------------|-------|
| LN-72045512 | Replacement Filters, Pack of 5 | |



Battery and Filter not included in warranty



Soft-Coated Triggers
with direction indicators

Ergonomically Balanced
in the hand with a loaded pipette

Universal Grip
for right- or left- handed users



Adapter
for small pipettes

Easy Identification

5 different colored nose cones help differentiate department, task or user; as well as a labeling area on the side of the unit



Magnetic & Screw Wall-Mount



3 Modes

Bright LEDs indicate 3 operation modes: high, low, & gravity



RF3000®

Packed with features for comfortable and accurate pipetting

The powerful pump offers rapid, yet sensitive aspiration and dispensing, has adjustable valves and is microprocessor controlled for constant reliability and accuracy

- Long-life li-ion battery, with quick charge—up to 20 hours use
- Toggle operation by pressing the large mode button
- An adapter can be fitted to the nose cone to stabilize small volume and thinner pipettes
- New USB connector ensures the appropriate charger is used with the unit
- Low-battery indicator
- Can be recharged when in use



Nose cone assembly only



Storage Stand
doubles as a charging station

Specifications

| | |
|-----------------------|--|
| Power | Micro-USB connector Rechargeable li-ion battery (replaceable) |
| Filter | PTFE/25mm dia/0.45 µm filter |
| Pipettes | 1–100 mL, glass and plastic |
| Pipette Nozzle | Polypropylene |



Package Includes

- RF3000® Pipette Controller
- Charging Stand and Wall-Mountable/Magnetic Holder
- Nozzle Adapter for use with Small Pipettes
- Li-Ion Rechargeable Battery
- Five Nose Cones: Black, Yellow, Red, Grey and Blue
- 5VDC Micro USB Plug
- 6 ft (183 cm) Cable
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



| Item No. | Volume Range mL | L x W x H in | L x W x H cm |
|-------------|--------------------|-----------------|--------------------|
| LN-72045502 | 1–100 | 9 x 8.5 x 5 | 22.9 x 21.6 x 12.7 |

Replacement parts

| | | | |
|-------------|---|----------|--|
| LN-72045512 | RF3000® Replacement Filters, Pack of 5 | | |
| LN-72045158 | RF3000® Universal Power Adaptor 5VDC Micro USB Plug | | |
| LN-72045517 | RF3000® Nose Cones Assorted* ■ ■ ■ ■ ■ | Set of 5 | |
| LN-72045750 | RF3000® Replacement Filter Adapter | | |
| LN-72045504 | RF3000® Replacement Pipette Silicone Adapter | | |

*Assorted colors (Black, Blue, Grey, Red and Yellow)



eCount

Rapid, accurate colony counter

Combines a multi-function electronic counter with a Sharpie® Pen for marking to prevent missing or double-counting colonies

- Light, ergonomic body is balanced for easy manipulation
- Select between an audible beep or an LED flash to acknowledge count
- Store up to 32 separate counts and display the total of all counts
- Count up/down verification
- Quick and easy to change between different colored markers to check different colony counts in the same view

Package Includes:

- eCount Colony Counter
- Black Sharpie® Pen (Page 16)
- Long-life button battery included
- Hand tally Adapter
- Weighted Stand



US Patent No. D564,385
EP Registered Design No. 705751-0001

| Item No. | L x W in | L x W cm |
|-------------|----------------|-------------|
| LN-72045163 | 7.4 x 1.2* | 17 x 3* |
| LN-72045161 | CR1225 Battery | |

*Measurements shown are with Sharpie® Pen inserted



Weighted Stand
for bench top use



Battery and Pen not
included in warranty



Lock Counter Off

to allow use of the marker without
registering more counts

Hand Tally Adapter

allows you to use the
eCount as a digital tally
counter





Mini-Light Box and Mini-Magnifier

Perfect accessories for the eCount or use as stand-alone items

Mini-Light Box



- LED Light source. Color temperature 5,300K +/- 270K
- Viewing area (125 x 100 mm) is perfect for quick screening of 96-well plates, petri dishes, electrophoresis gels, slides, and transparencies
- Powered by six 1.5V AAA batteries (included), not included in warranty

Mini-Magnifier

- Swiveling molded lens 1.75x magnification
- Enlarges objects for easier viewing

| Item No. | Description | L x W x H in | L x W x H cm |
|-------------|-----------------|-----------------|--------------------|
| LN-72045154 | Mini-Light Box* | 6.2 x 5.9 x 0.6 | 15.8 x 15 x 1.4 |
| LN-72045146 | Mini-Magnifier | 6.3 x 5.3 x 8.1 | 16.1 x 13.5 x 20.7 |

*Batteries included



Hand Tally Counter

Use for counting blood cells, bacterial colonies, drops of liquid or any repetitive event

- Tallies up to 9999
- Features a quick-reset knob
- Lens for easy reading, and finger ring for a secure grip

| Item No. | Dia. x W in | Dia. x W cm |
|-------------|----------------|----------------|
| LN-72042723 | 1.8 x 1.6 | 4.6 x 4.1 |



Benchtop Counters

Mechanical, benchtop counters are perfect for blood cells (label included), bacterial colonies or repetitive counting

- Each key records up to 999
- Blood cell key labels
- Total window keeps track of total number of strokes made on all other keys; bell rings when total reaches 100

| Item No. | No. of Keys | L x W x H in | L x W x H cm |
|-------------|-------------|-----------------|-----------------|
| LN-72042717 | 6-Key | 9.8 x 3 x 2.2 | 25 x 7.6 x 5.6 |
| LN-72042720 | 9-Key | 12.5 x 3 x 2.2 | 32 x 7.6 x 5.6 |

Laboratory Notebooks

Keep information safe and well organized with Laboratory Notebooks

Record research data, experimental procedures, design and development ideas and more

- Quality Smyth Sewn binding holds pages securely and lies flat when opened
- Made of acid-free paper, pages are numbered and include space for title, date, project number, book number, date and witness signature
- User data, notebook instructions, table of contents, periodic table, and conversion tables are included in every book

Keep Information Safe
with a durable hardbound, water-resistant cover




Easy Identification
with a foil stamped title

Marker Ribbon
to keep your place

| Item No. | Color | Pages | Page Format. | L x W in | L x W cm |
|-------------|-------|-------|--------------|-------------|-------------|
| LN-72042748 | Black | 100 | Grid | 8.6 x 11.1 | 21.8 x 28.1 |
| LN-72042750 | Blue | 200 | | | |
| LN-72042753 | Red | 200 | Lined | | |
| LN-72042746 | Black | 100 | | | |

Case of 12



2 Page Formats

- Grid: 6 x 6 mm grid size
- Lined: college-ruled pages for easy report writing



Twin Tip & Fine Tip Sharpie® Permanent Pens

Ink dries quickly and resists both fading and water

- Pen has a fine and an ultra fine tip and comes in a bright selection of colors
- Ink dries quickly and resists both fading and water
- AP certified

| Item No. | Description | Color |
|-------------|--------------------------------|----------|
| LN-72041592 | Twin Tip Marker Set, Pack of 8 | Assorted |
| LN-72041589 | Fine Point Marker, Pack of 12 | Black |

*Assorted colors (Black, Turquoise, Navy, Lime, Orange, Purple, Sky Blue, Berry)
Sharpie® is a registered trademark of Sanford, A Newell Rubbermaid Company

Lab Alert® Clocks & Timers



Pocket Timer/ Stopwatch/Clock

- Four independent timing channels with visual and audible alarms
- Use as a stopwatch, a count up/down timer, or a clock
- Includes spring clip and magnetic back



Triple Timer/Clock

- Three channels to count up/down with
- Simultaneous display, time-out, and 60 second audible alarm
- Channels may be started independently or simultaneously
- Timer automatically counts up after counting down to zero
- Memory recall



4-way Timer with Whiteboard/Pen

- Clock feature count up or down to 24 hours 59 minutes
- 4 timer channels
- Stands up or attaches to a magnetic surface (magnets on back of unit)
- Whiteboard and marking pen included
- Sounds when time expires

| Item No. | No. of Channels | Timing Capacity | Count Direction | Display Digits, H mm | Battery Included | L x W x H in | L x W x H cm |
|---------------|-----------------|-----------------|-----------------|----------------------|------------------|-----------------|------------------|
| a LN-72042290 | Four | 100 hr. | Up/Down | 6 | One LR44 1.5V | 2.4 x 2.8 x 1 | 6.1 x 7.1 x 2.5 |
| b LN-72042291 | Three | 20 hr. | Up/Down | 10.5 | One AAA | 3.4 x 2.6 x 0.9 | 8.6 x 6.6 x 2.2 |
| c LN-72042292 | Four | 24 hr. | Up/Down | 10 | Two AAA | 4.8 x 3.4 x 0.9 | 12.3 x 8.6 x 2.1 |
| d LN-72042294 | One | 24 hr. | Up | 8.5 | One LR44 1.5V | 3.4 x 3.1 x 0.8 | 8.7 x 7.8 x 2.1 |
| e LN-72042303 | Four | 99 hr. | Down | 14 | Two LR44 1.5V | 2.9 x 2.3 x 0.5 | 7.3 x 5.9 x 1.3 |



Multi-function Stopwatch/Clock

- Quick start, quick stop timer features
- time, split, hourly chime, calendar, alarm
- 100 cm neck cord



Four Channel Timer/Clock

- Independently or simultaneously set time on four channels
- Audible and visual alarms signify elapsed time
- Easily toggle between channels and clock with top buttons
- Includes spring clip with magnetic back

Lab Drying Rack

Bringing color to a mundane task

Most Versatile Rack in the market


- Channels designed for more useful configurations
- Large and Small pegs for all glassware varieties
- One of a kind Drying Basket to collect smaller items

Designed to protect against breakage

- Peg angle allows for better hold - especially large, heavy pieces
- Patent Pending SureGrip Hook securely grips glassware into place

Included:

- 72 Pegs ~ 36 Small (Green) ~ 36 Large (Blue)
- 2 Drying Baskets
- Drainage Tubing
- Mounting Hardware



Drying Basket

to keep small items and slotted to hold pipets safe while drying



SureGrip Hook

Patent-pending design securely grips glassware into place




Use on counter or wall

Stands alone on integrated hinged stand. Mounts easily on wall with included guide and hardware

Can hold up to 77 lbs
(34.9 kg)

Rack made of durable polystyrene



Lab Drying Rack

U.S. Patent No. D923275, D936925, D951573, D951574

| Item No. | L x W x H in | L x W x H cm |
|-------------|-----------------|------------------|
| LN-72042845 | 17.4 x 26 x 5.5 | 44.1 x 66 x 13.9 |

Replacement Parts

| Item No. | Description |
|-------------|--|
| LN-72042849 | Peg 10 Piece Kit, (5) Small, (5) Large |
| LN-72042847 | Drying Basket Kit, 2 ea |
| LN-72042851 | Replacement Draining Set, Tube and Elbow |

Workstation Lab Mat

Safely facilitate workflow processes

- Protect counter tops from hot items up to 200°C (392°F)
- Lip around the edge can hold approx. 250 mL of liquid
- Small footprint, 12 x 14.9 x 0.1 in (30.5 x 37.9 x 0.2 cm)



Gridded Surface

for tube staging, preventing movement of objects

Reversible

Side One: combination gridded and smooth surface

Side Two: smooth surface only



| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--|-----------------|-------------------|
| LN-72045309 | Blue ■ | 12 x 14.9 x 0.1 | 30.5 x 37.9 x 0.2 |
| LN-72045310 | Purple ■ | | |

Silicone Lab Mat

Large working space

The lab mat can be used with items such as centrifuges, vortex mixers, hot plates, stirrers, pipettes and so much more

- Mat features designated magnetic area to keep stir bars in place
- Protect countertops from hot items up to 200°C (392°F)




Raised Edges

enable liquid retention of approximately 300 mL of fluid



Large Working Surface

accommodates benchtop equipment, weigh boats, stir bars, scales, etc.



| Item No. | Side #1 Color | Side #2 Color | L x W x H in | L x W x H cm |
|-------------|--|---|-----------------|-----------------|
| LN-72041880 | Yellow ■ | Grey/Blue ■ ■ | 14 x 23.5 x 0.1 | 35 x 60 x 0.2 |
| LN-72041881 | Purple ■ | Purple/Grey ■ ■ | | |
| LN-72041882 | Blue ■ | White/Blue ■ ■ | | |

Blood Collection Tray

Light weight, heavy on organization

Organize your phlebotomy, I.V. start or catheter insertion supplies with this spacious, robust and easy to clean carry tray

- Lightweight easy-to-clean ABS construction
- Customize the compartments to suit your requirements with the dividers provided
- Comfort grip on durable acetal plastic handle tested to hold up to 10 kg

| Item No. | L x W x H in | L x W x H cm | Handle Height in | Handle Height cm |
|-------------|-------------------|-------------------|---------------------|---------------------|
| LN-72041454 | 16.4 x 12.2 x 3.3 | 41.9 x 31.1 x 8.5 | 9.4 | 24.1 |

Functional

Comes with a 17 mm well,
60-place, tube rack



Totemax® Blood Collection Tray

Keep off the floor, work top or bed

Unique phlebotomy tray design and patented hook feature allows the tray to be on or over most bed rails, foot-boards and handrails found in a typical patient room

- Made from food grade HDPE resins for a durable non-breakable design
- Deep tray design helps reduce spills



Hooks Lock

in place when not in use

Hangs on bed rails,
footboards and handrails



| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--------|-----------------|--------------------|
| LN-72041136 | Blue ■ | 17 x 11 x 11 | 43.2 x 27.9 x 27.9 |



Droplet Blood Collection Tray

Designed with phlebotomists in mind

Keeping a clean and organized tray has never been easier. Wash the ABS tray as needed, and replace the insert periodically for a fresh start

- Multiple compartments keep a variety of supplies neatly at hand, and organized
- Built in side pockets offer additional storage space
- Handle includes soft rubber grip to reduce hand and finger strain
- Each tray includes either a 13 or 16 mm Heathrow OneRack® tube rack (page 34) and two disposable polystyrene inserts

Ratcheting Handle

folds out of the way for easy access, even over tall sharps containers



Supply Drawer

Side Pockets

Tray Inserts

come in two configurations and include an extra-deep compartment for a sharps container



Disposable Insert A
for 13 mm tube rack



Disposable Insert B
for 16 mm tube rack

Droplet Configurations

US Patent No. D544703

| Item No. | Tray Insert | Tube Rack mm | L x W x H in | L x W x H cm | Handle Height | |
|-------------|-------------|-----------------|-------------------|--------------------|---------------|----|
| | | | | | in | cm |
| LN-72041973 | HS2201A | 13 | 16.8 x 10.4 x 6.7 | 42.9 x 26.5 x 17.2 | 11 | 28 |
| LN-72041976 | HS2201B | 16 | | | | |

Droplet Inserts

Pack of 4

| Item No. | Description | L x W in | L x W cm |
|-------------|---------------------------------|-------------|-------------|
| LN-72041979 | Tray insert for 13 mm tube rack | 14.2 x 8.2 | 36 x 21 |
| LN-72041982 | Tray insert for 16 mm tube rack | | |



Keeps contents clean and secure

Droplet Lid

Cover for the Droplet blood collection tray

Protect the contents of the Droplet Tray with this Droplet lid. Custom designed sturdy lid is a snug fit, protects contents and provides access to the tray handle in its vertical position

- Made of polystyrene
- Can be secured with cable ties for extra security
- Stackable



| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--------------------------------|-----------------|-----------------|
| LN-72041985 | White <input type="checkbox"/> | 16.5 x 9.4 x 8 | 42 x 24 x 20.5 |

Carry handle folds

neatly into body space allowing the units to be stacked

Molded Slots

to hold tubes while working

**Silicone Water-tight Seal**

and three lid clasps provide secure closure

**Lid Opens a Full 180°**

to expedite loading and emptying

DuraPorter® Specimen or Sample Transport Box

Ideal for transporting biological and clinical samples, as well as instruments and products that you want to keep clean and dry under testing environmental conditions

The Duraporter® transparent design is OSHA 29 CFR 1910.1031 and UN3373 compliant

- Polycarbonate construction and available in three colors for user/owner identification
- Use the separators provided to create your own compartment sizes
- Accommodates 13 mm, 16 mm and 20 mm tubes with our Half and Full Size OneRack®



| Item No. | Color | | L x W x H in | L x W x H cm |
|-------------|---------------------------|---|-----------------|-----------------|
| LN-72041193 | Clear with Blue Handle | ■ | | |
| LN-72041229 | Red with Red Handle | ■ | 15 x 7.7 x 6.3 | 38 x 19.6 x 16 |
| LN-72041232 | Yellow with Yellow Handle | ■ | | |

DuraPorter® XL

High-Capacity Sample Transport Tote

Your reliable solution for the safe and convenient transportation of laboratory samples, instruments, and valuable items. Crafted to meet the stringent demands of clinical and research laboratories, this oversized transport container ensures the protection and preservation of your cargo

Built to Last and Designed with the User in Mind

- Sturdy and spacious, ideal for transporting a large volume of samples
- Transparent design helps comply with OSHA 29 CFR 1910.1030 standards without having to re-label containers and enabling the user to visibly check sample integrity and meets UN3373 secondary packaging transport requirements
- Equipped with articulated handles featuring silicone grips for comfortable and secure handling

Emphasizing Security and Safety


- A silicone seal ensures a tight enclosure, safeguarding samples from most environmental elements
- Dual latches with security tag feature guarantee a secure closure, providing tamper-evidence and ensuring the integrity and authenticity of your samples

Soft Grip Handle
for comfortable handling

Removable Lid
for ease of loading and unloading

4-Point Handle Security
ensures durable transport

Translucent Design
allows for easy content identification



Silicone Seal
for secure content enclosure

Dual Latches
with tag out feature for added security



DuraPorter® XL

Includes One 120962 MultiTube Rack Insert

| Item No. | Color | Dimensions (not including handle) in | Dimensions (not including handle) cm | UOM |
|-------------|--------|--------------------------------------|--------------------------------------|------|
| LN-72041235 | Blue | 16 x 16.25 x 7 | 40.64 x 41.28 x 17.78 | 1 ea |
| LN-72041236 | Red | 16 x 16.25 x 7 | 40.64 x 41.28 x 17.78 | 1 ea |
| LN-72041237 | Yellow | 16 x 16.25 x 7 | 40.64 x 41.28 x 17.78 | 1 ea |


DuraPorter® XL

Tube Rack Inserts Accessories



Duraporter® XL Large Tube Insert

Holds 100 tubes with diameter up to 18 mm
Most 10 mL blood tubes & 15 mL conical tubes
Array 10 x 10
Pack of 3



Duraporter® XL Small Tube Insert

Holds 168 tubes with diameter up to 13 mm
Most 3.5 - 7.5 mL blood tubes
Array 12 x 14
Pack of 3



Duraporter® XL MultiTube Rack Insert

Holds tube with various diameters: 13, 18 & 30 mm
Most 15 mL & 50 mL conical tubes
and 3.5 mL - 10 mL blood tubes
Pack of 3

DuraPorter® XL Tube Rack Inserts

| Item No. | Tubes | Cut-out Diameter | UOM |
|-------------|---|--|-----------|
| LN-72041238 | 15 mL, Large blood tube | 18 mm | Pack of 3 |
| LN-72041239 | Small blood tubes | 13 mm | Pack of 3 |
| LN-72041240 | 50 mL, 15 mL, Small & Large blood tubes | 30 mm (center), 18 mm (large), 13 mm (small) | Pack of 3 |

Aluminum Weighing Boats & Dishes

Dishes for general-purpose weighing, dispensing, storage, food processing, evaporating or drying



AC-74160070



AC-74200070



AC-74300082



AC-74160028

Aluminum Weighing Boats

| Item No. | Volume mL | Dimensions in | Dimensions mm | Packaging | Tab | Wall Finish |
|-------------|--------------|------------------|------------------|-----------|-----|-------------|
| AC-74200140 | 140 mL | 4.1 x 0.8 | 103 x 20 | 100 | | Smooth |
| AC-74200070 | 80 mL | 2.9 x 0.7 | 73 x 18 | 100 | | Smooth |
| AC-74300087 | 70 mL | 2.8 x 0.7 | 71 x 17 | 100 | | Crinkle |
| AC-74160070 | 50 mL | 2.8 x 0.6 | 70 x 16 | 100 | x | Crinkle |
| AC-74200057 | 42 mL | 2.2 x 0.6 | 57 x 16 | 100 | | Smooth |
| AC-74160060 | 39 mL | 2.2 x 0.6 | 57 x 14 | 100 | x | Crinkle |
| AC-74300082 | 20 mL | 2 x 0.5 | 52 x 13 | 100 | | Crinkle |
| AC-74160045 | 20 mL | 1.7 x 0.5 | 43 x 12 | 100 | x | Crinkle |
| AC-74160028 | 5 mL | 1.1 x 0.4 | 28 x 10 | 100 | x | Crinkle |

Aluminum Weighing Boats

- Sturdy yet lightweight boats feature a smooth, flat bottom for sample handling

Aluminum Weighing Dish

- Designed for moisture determination balances with slightly tapered walls for stacking



AC-74120125



AC-74120095

Aluminum Weighing Dish

| Item No. | Volume mL | Dimensions in | Dimensions mm | Packaging | Tab |
|-------------|--------------|------------------|------------------|-----------|-----|
| AC-74120095 | 45 mL | 3.5 x 0.3 | 90 x 7 | 50 | |
| AC-74120125 | 40 mL | 4 x 0.3 | 101 x 8 | 50 | |
| AC-74120065 | 15 mL | 2.7 x 0.2 | 70 x 6 | 50 | |
| LN-72041516 | 1 mL | 0.98 x 0.3 | 25 x 7 | 50 | |
| 120933FR | 1 mL | 0.98 x 0.3 | 25 x 7 | 50 | |
| LN-72041517 | 0.25 mL | 0.5 x 0.1 | 13 x 3.5 | 100 | x |
| 120940FR | 0.25 mL | 0.5 x 0.1 | 13 x 3.5 | 100 | x |



AC-74120065



LN-72041517

LevGo® smartPans®

Ideal for microbalances, organic micro elementary and thermal analysis, and pharmaceutical use

- Made of Ultra High Vacuum (UVH) Foil - "A" wettable and completely free from oils normally found on other grades of aluminum foil
- Foil meets ASTM B-479, specification 3.1.4 and 10.3.1
- "Cahn style" accessory pans are compatible with most commercially available analyzers and microbalances



| Item No. | Dia. x H in | Dia. x H mm | UOM |
|-------------|----------------|----------------|-------|
| LN-72041501 | 1 x 0.1 | 25 x 3.6 | 20/Pk |

Weighing Papers

Smooth surface for good sample recovery

- Nitrogen-free
- Strong enough to support dense material



Pack of 500

| Item No. | L x H in | L x H cm |
|-------------|-------------|-------------|
| AC-43760707 | 3 x 3 | 7.6 x 7.6 |
| AC-43761010 | 4 x 4 | 10.2 x 10.2 |
| AC-43761515 | 6 x 6 | 15.2 x 15.2 |



Functional in the Lab

Available in a small and medium in both brown for light colored samples and white for dark samples



The mark of responsible forestry

LevGo® smartBoats®

A truly green solution to sample weighing

Made from FSC® Certified (FSC® Mix Credit) paper (Eco-IN) and when disposed of in the standard waste stream they will naturally biodegrade (Eco-OUT)

- A cost-effective solution to reduce plastics in the lab
- Made of grease-proof paper that meets FDA 21 CFR 176.170 (components of paper in contact with aqueous and fatty foods) and with 21 CFR 176.180 (contact with dry foods)
- Control contents like a polystyrene weigh boat because they are easier to pour
- An alternative to regular weighing papers because they are easier to handle

| Item No. | Color | Size | Packaging | L x W x H in | L x W x H cm |
|-------------|---|--------|-------------|-----------------|-----------------|
| AC-72630040 | Brown <input checked="" type="checkbox"/> | Small | Pack of 500 | 1.4 x 1.4 x 0.6 | 3.5 x 3.5 x 16 |
| AC-72650040 | White <input type="checkbox"/> | | | | |
| AC-72630050 | Brown <input checked="" type="checkbox"/> | Medium | Pack of 500 | 1.9 x 1.9 x 0.6 | 4.8 x 4.8 x 14 |
| AC-72650050 | White <input type="checkbox"/> | | | | |
| AC-72630090 | Brown <input checked="" type="checkbox"/> | Large | Pack of 250 | 3.8 x 3.8 x 1 | 9.5 x 9.5 x 2.5 |
| AC-72650090 | White <input type="checkbox"/> | | | | |

smartBoats™ is a trademark of LevGo Inc.

LevGo® smartSpatula®

Designed to fit perfectly with standard
0.2, 0.5, 1.5/2.0, 15 and 50 mL tubes

- Made of virgin, food-grade polypropylene
- Ideal for dispensing or sample collection of powders or soft solids
- Will not shatter when used with liquid nitrogen frozen samples
- Hollow core minimizes heat transfer—good for hot and cold liquids
- Antistatic, sterile and Eco-Friendly versions available
- Complies with FDA standards for food contact



a



b



c



d



e

AUTOCLEAVABLE
121° C

STERILE EO

Always there
and always
clean!

Lab tools that scoop, stir and scrape.



| Item No. | Description | Color | Suggested Tube mL | Packaging | L x W in | L x W cm | Spoon (Flare) in | Spoon (Flare) mm | Scoop/Pick L x W in | Scoop/Pick L x W mm | Taper |
|---------------|-------------------|---|-------------------|---------------|-------------|-------------|---------------------|---------------------|------------------------|------------------------|-------|
| a LN-72041103 | Micro Antistatic | Opaque <input type="checkbox"/> | 0.2 and 0.5 | 300/Bulk | 5.5 x .2 | 14 x .4 | — | — | 4.3 | 11 (V-Tip) | 15° |
| b LN-72041112 | Eco-Friendly | Opaque <input type="checkbox"/> | 15 and 1.5 | 300/Bulk | 5.5 x .3 | 14 x .7 | 9.8 x 5.5 | 25 x 14 | 11.8 x 1.2 | 30 x 3 | |
| c LN-72041097 | Standard | Blue <input checked="" type="checkbox"/> | | 300/Bulk | | | | | | | |
| LN-72041100 | Standard | Opaque <input type="checkbox"/> | | 300/Bulk | | | | | | | |
| d LN-72041109 | Std. EtO Sterile | Opaque <input type="checkbox"/> | 15 and 1.5 | 100/10 x 10pk | 8.3 x .3 | 21 x .7 | 9.8 x 5.5 | 25 x 14 | 11.8 x 1.2 | 30 x 3 | |
| LN-72041110 | Std. EtO Sterile* | Opaque <input type="checkbox"/> | | 100/10 x 10pk | | | | | | | |
| e LN-72041106 | Macro | Green <input checked="" type="checkbox"/> | 50 | 150/Bulk | 12.2 x .4 | 31 x 1 | 16.5 x 9.3 | 42 x 23.5 | 11.8 | 30 | 35° |

*Non Toxic/Non Pyrogenic smartSpatula™ is a trademark of LevGo Inc.

Weighing Boats

Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss

- Select white, black, or blue for best color contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Antistatic versions available making the boats nonreactive with most substances and suitable for weighing static-affected samples
- Wide flat bottoms ensure stability and prevent tipping
- Rounded corners allow simple flexing for controlled pouring of liquid, powdered or granular samples
- Can be used with microwave moisture analyzers
- Made from FDA food grade polystyrene



Standard Weighing Boat

Weigh Boats

Pack of 500

| Shape | Antistatic Material Item No. | Color | Size | Volume mL | L x W x H in | L x W x H cm | |
|-----------|------------------------------|---|---|-----------|-----------------|-------------------|-----------------|
| Standard | AC-70240040 | White <input type="checkbox"/> | Small | 10 | 1.8 x 1.8 x .3 | 4.6 x 4.6 x .8 | |
| | AC-70220040 | Blue <input checked="" type="checkbox"/> | | | | | |
| | AC-70260040 | Black <input checked="" type="checkbox"/> | | | | | |
| | Standard | AC-70240080 | White <input type="checkbox"/> | Medium | 85 | 3.4 x 3.4 x .9 | 8.5 x 8.5 x 2.4 |
| | | AC-70220080 | Blue <input checked="" type="checkbox"/> | | | | |
| | | AC-70260080 | Black <input checked="" type="checkbox"/> | | | | |
| | | AC-70240140 | White <input type="checkbox"/> | | | | |
| | | AC-70220140 | Blue <input checked="" type="checkbox"/> | | | | |
| | AC-70260140 | Black <input checked="" type="checkbox"/> | | | | | |
| Hexagonal | AC-70540025 | White <input type="checkbox"/> | Small | 6 | 1.5 x 1 x .4 | 3.8 x 2.5 x .9 | |
| | AC-70540052 | | Medium | 46 | 2.5 x 2 x .7 | 6.4 x 5.1 x 1.9 | |
| | AC-70540085 | | Large | 194 | 4 x 3 x .9 | 10.2 x 7.6 x 2.2 | |
| | AC-70540093 | | X-Large | 364 | 5 x 3.5 x 1.2 | 12.7 x 8.9 x 3.1 | |
| Pour-boat | AC-70440044 | White <input type="checkbox"/> | Small | 14 | 1.7 x 2.3 x .5 | 4.3 x 5.8 x 1.3 | |
| | AC-70440090 | | Medium | 137 | 3.3 x 5.2 x 1 | 8.3 x 13.2 x 2.6 | |
| | AC-70440120 | | Large | 259 | 4.3 x 7.2 x 1 | 10.8 x 18.3 x 2.6 | |
| Square | — | White <input type="checkbox"/> | Medium | 51 | 3.5 x 3.5 x .3 | 9 x 9 x .7 | |
| Diamond | AC-70340055 | White <input type="checkbox"/> | Small | 5 | 2.1 x 1.3 x .2 | 5.3 x 3.3 x .6 | |
| | AC-70360055 | Black <input checked="" type="checkbox"/> | | | | | |
| | BR-70155548 | White <input type="checkbox"/> | Large | 30 | 3.1 x 2.2 x .6 | 7.8 x 5.6 x 1.4 | |
| | AC-70360078 | Black <input checked="" type="checkbox"/> | | | | | |
| | BR-70155550 | White <input type="checkbox"/> | X-Large | 100 | 4.5 x 3.3 x 0.7 | 115 x 85 x 17 | |
| | AC-70360120 | Black <input checked="" type="checkbox"/> | | | | | |



Standard Weighing Boat



Square Weighing Boats



Hexagonal Weighing Boats



Diamond Weighing Boats



Pour-boat Weighing Dishes

32 Place PCR Rack

Designed for pre- and post-PCR sampling, streamline your PCR workflow with the 32-Place PCR Rack. This versatile rack accommodates 0.2 mL tubes both individually and in strips of 8, as well as 32-place assay plates (8x4)

Space Saving Design

- Small and compact design is stackable, saving valuable space in the freezer
- Right sized when processing smaller batches
- Lid hinge design allows lid to sit at 180 degrees, clear of sample processing, or can be removed entirely

Reusable Convenience

- Transparent lid protects samples and allows for quick assessment of contents
- Embossed numbers and letters on the wells provide easy identification and tracking of samples
- Light-weighted construction made durable for reusability

Pack of 10

| Item No. | Tubes mL | No of Wells | L x W x H in | L x W x H cm |
|-------------|-------------|-------------|-----------------|------------------|
| LN-72042134 | 0.2 | 32 | 3.8 x 6.9 x 1.2 | 9.8 x 17.5 x 3.1 |

Made from 100% recyclable Polypropylene

Stores down to -80°C



Hinged Lid
stays clear of sample processing

Embossed Letters
& Numbering



PCR Tube Racks 96-Well

Perfect for pre- or post-PCR sampling



- Made of polypropylene
- Accommodates individual tubes or strips
- Removable hinged lid lies flat
- Alphanumeric identification for quick location and retrieval of tubes

Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|-----------|-------------|-------|
| LN-72042153 | Assorted* | LN-72042159 | Blue |
| LN-72042156 | Natural | LN-72042162 | Green |

*Assorted colors (Blue, Green, Pink, Yellow, Orange)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|------------------|
| 0.2 | 96 | 5.2 | 5 x 3.9 x 1.2 | 12.7 x 9.9 x 3.0 |



1.5 mL Tubes
12 wells on each side

Reversible
for 0.2 mL tubes
or for 0.5 mL tubes

PCR Tube Racks Reversible

Removable hinged lid fits both sides



- Made of polypropylene
- PCR side of the rack holds 8 or 12 x 0.2 mL tube strips or individual tubes
- Reverse side holds 0.5 mL tubes
- Both sides of the rack have 12 wells for 1.5 mL tubes

Pack of 5

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in (cm) |
|-------------|-------|-------------|-------|-----------------|----------------------|
| LN-72042078 | | 0.2 | 168 | 5.5 | 8.1 x 5.2 x 2.1 |
| | | 0.5 | 40 | 7.6 | (20.6 x 13.1 x 5.3) |
| | | 1.5/2.0 | 12 | 11.4 | |

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



Alphanumeric Identification
for quick location and retrieval of tubes

Microtube Storage Boxes

- Hold 1.5 and 2.0 mL microcentrifuge tubes and cryogenic vials




Stylish & Identifiable

- Made of translucent, vibrantly colored polypropylene allows for visibility inside box and color coding identification
- Tube capacity is indicated on the outside of the box within a distinctive wave design
- Quick visual identification with uniform wells and alphanumeric location identification along with a frosted labeling area

Robust

- Easy open and secure lid fitted with either friction fit or snap clasp designs
- Stackable, saving valuable space in the lab



50-Well Microtube Storage Boxes

Pack of 5

| Item No. | Color | |
|-------------|-----------|--|
| LN-72041154 | Natural | <input type="checkbox"/> |
| LN-72041157 | Assorted* | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> |

**Assorted colors (Blue, Green, Purple, Yellow, Orange)*

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|------------------|
| 1.5/2.0 | 50 | 12.3 | 5.6 x 3.6 x 2.2 | 14.2 x 9.2 x 5.6 |



81-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|----------------------------------|-------------|--|
| LN-72041160 | Natural <input type="checkbox"/> | LN-72041172 | Yellow <input type="checkbox"/> |
| LN-72041163 | Blue <input type="checkbox"/> | LN-72041175 | Orange <input type="checkbox"/> |
| LN-72041166 | Green <input type="checkbox"/> | LN-72041178 | Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> |
| LN-72041169 | Pink <input type="checkbox"/> | LN-72041181 | Black <input type="checkbox"/> |

**Assorted colors (Blue, Green, Yellow, Orange, Pink)*

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|-----------------|
| 1.5/2.0 | 81 | 12.3 | 5.1 x 5.1 x 1.9 | 13 x 13 x 4.7 |



100-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|--|----------|---------------------------------|
| LN-72041184 | Natural <input type="checkbox"/> | 120202 | Green <input type="checkbox"/> |
| LN-72041187 | Black <input type="checkbox"/> | 120203 | Yellow <input type="checkbox"/> |
| LN-72041190 | Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> | 120204 | Orange <input type="checkbox"/> |
| LN-72041346 | Blue <input type="checkbox"/> | | |


**Assorted colors (Purple, Blue, Green, Yellow, Orange, Pink)*

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|-------------------|
| 1.5/2.0 | 100 | 12.3 | 5.6 x 6 x 2.2 | 14.1 x 15.1 x 5.7 |



Tabs on each end
of rack snap lid into place

Raised Rings
around the 1.5 to 2.0
mL wells make tubes
easier to remove



Molded Grid Lines

0.5 mL wells separated by molded-in
grid lines, making them easier to see



Reversible Tube Racks 96-Well

Easy-to-grasp tabs for sample handling and storage

96 wells for 0.5 mL tubes on one side and 96 wells for 1.5/2.0 mL tubes on the reverse side

Optimum Function & Identification

- Wells arranged in standard 8 x 12 array
- Molded in-grid lines separate 0.5 mL wells making them easier to see
- Raised rings around the 1.5 mL/2.0 mL wells make tubes easier to remove
- Alphanumeric identification for quick location and retrieval of tubes

Secure lid

- Translucent lid snaps onto rack tabs for secure, visible storage
- Spacious enough to fit most screw-top micro-centrifuge tubes
- Stackable, saving valuable space



Pack of 5


| Item No. | Color | Item No. | Color |
|-------------|-----------|-------------|--------|
| LN-72042132 | Assorted* | LN-72042141 | Green |
| LN-72042135 | Natural | LN-72042144 | Pink |
| LN-72042138 | Blue | LN-72042147 | Yellow |

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|-----------------|
| 1.5/2.0 | 96 | 11 | 9.6 x 4.7 x 2 | 24.4 x 12 x 5.1 |
| 0.5 | | 7.6 | | |

Easy Identification

Alphanumeric well identification



Rota-Rack® Duo Tube Rack

Dilution and concentration studies made easy

Transfer solutions between 0.2 mL PCR tubes and 0.5 and 1.5 mL microcentrifuge tubes

- Purple section holds 48 x 0.5 mL and 48 x 1.5/2.0 mL microcentrifuge tubes on the reverse side
- Blue section holds 96 x 0.2 mL tubes and 48 x 1.5/2.0 mL tubes on the reverse side



| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-----------------|-------------|-------|-----------------|------------------|-------------------|
| LN-72042642 | Blue/ Purple | 0.2 | 96 | 5 | 11.6 x 4.6 x 1.3 | 29.6 x 11.8 x 3.4 |
| | | 0.5 | 48 | 8 | | |
| | | 1.5/2.0 | 96 | 11 | | |

Two Tier Microtube Rack 20-Well

Interlock multiple racks to increase tube storage

Improved Experience

- Tiered rack allows for improved staging and visibility for easier, more accurate sample filling and retrieval
- Made of easy to clean polypropylene
- Easily remove tubes using finger slot on front of the rack
- Embossed Index Numbers for quick sample ID and retrieval



Pack of 6

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-------------|-------|-----------------|-----------------|-----------------|
| LN-72042675 | Blue | 1.0 | 20 | 10 | 8.3 x 3.5 x 2 | 21.1 x 9 x 5 |
| | | 1.5 | | | | |
| | | 2.0 | | | | |



Microtube Racks 80-Well

Popular fraction-collector-plate style rack

Compact 5x16 array makes filling and transfer easy

- Frosted labeling area on bottom of rack
- Embossed index numbers for quick sample ID and retrieval
- Made of polypropylene



Pack of 5

| Item No. | Color | Part No. | Color |
|-------------|---------|-------------|-----------|
| LN-72042612 | Natural | LN-72042624 | Yellow |
| LN-72042615 | Blue | LN-72042627 | Orange |
| LN-72042618 | Green | LN-72042630 | Assorted* |
| LN-72042621 | Pink | LN-72042633 | Red |

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|----------|-------|--------------|---------------|------------------|
| 1.5/2.0 | 80 | 11 | 8.9 x 2.6 x 1 | 22.5 x 6.7 x 2.8 |



Storage Racks for 80-well Racks

Strong, clear acrylic construction

These racks help you keep large numbers of microtubes organized on your benchtop.

- Holds up to eight 80-well microtube racks filled with either 1.5 or 2.0 mL microtubes
- Plenty of room to slide the racks in and out

| Item No. | Description | Color | L x W x H in | L x W x H cm |
|-------------|-------------|-------|------------------|--------------------|
| LN-72042036 | Vertical | Clear | 6.2 x 9.3 x 9.2 | 15.7 x 23.6 x 23.4 |
| LN-72042039 | Horizontal | Clear | 12.2 x 9.3 x 4.8 | 31 x 23.6 x 12.2 |

Microtube Racks

Durable, Simple to Use

- Made from reinforced nylon offering superior chemical and heat-resistance over most other plastic based racks
- Handles provide easy and safe transportation
- Stable design can be used in water baths
- Alphanumeric identification for quick location and retrieval of tubes



| Item No. | Color | Tubes mL | Wells | Well Dia. mm | Packaging | L x W x H in | L x W x H cm |
|-------------|-------|----------|-------|--------------|-----------|-----------------|-----------------|
| LN-72041394 | Blue | 1.5/2.0 | 24 | 11.1 | 8 | 5.6 x 3.5 x 2 | 14.3 x 9 x 5.2 |
| LN-72041400 | | | 96 | | 4 | 8.3 x 5.2 x 2.1 | 21 x 13.1 x 5.4 |



Rota-Rack® Tube Racks

Get organized with this interlocked rack

Each module is interlocked and adjustable with just a twist allowing you to work with multiple sized tubes in the same rack

- Suitable for variety of culture, centrifuge, micro-centrifuge, round, flat or conical bottomed tubes
- Two versions to accommodate small and large tube sizes



| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|--------|-------------|-------|-----------------|-----------------|-----------------|
| LN-72042636 | Yellow | 50 | 2 | 29 | 8.8 x 4.3 x 4.3 | 22.4 x 11 x 11 |
| | Green | 10-15 | 8 | 20 | | |
| | Blue | 15-20 | 6 | 16 | | |
| | Purple | 10-13 | 10 | 13.5 | | |
| LN-72042639 | Yellow | 15 | 6 | 16.4 | 8.2 x 3.9 x 3.9 | 21 x 10 x 10 |
| | Green | 1.5/2.0 | 9 | 11.2 | | |
| | Blue | 0.5 | 12 | 8 | | |
| | Purple | 0.2 | 32 | 6.5 | | |



Rotate and Customize
to suit your application needs

Cube Tube Rack

Interlocking multi-tube racks

- Made of polypropylene
- Allows you to work with different sized tubes all in one unit



| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|--------|-------------|-------|-----------------|-----------------|--------------------|
| LN-72042645 | Blue | 50 | 4 | 30 | 4.3 x 4.3 x 4.5 | 10.8 x 10.8 x 11.3 |
| | Green | 15 | 10 | 17 | | |
| | Pink | 3-7 | 12 | 12.5 | | |
| | Yellow | 1.5/2.0 | 16 | 11.5 | | |

Pack of 5

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



4-Way Tube Rack

Link together to meet your application needs

Versatile

- Use as standalone unit or infinite units can be linked to create configurations to fit number of wells and tube sizes needed

Robust Functional Design

- Molded alphanumeric identification for quick location and retrieval of tubes
- Made from polypropylene that is chemical resistant and easy to clean
- Racks can be used in water baths and placed in freezer



Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|----------------------------------|-------------|------------------------------------|
| LN-72042591 | Natural <input type="checkbox"/> | LN-72042603 | Yellow <input type="checkbox"/> |
| LN-72042594 | Blue <input type="checkbox"/> | LN-72042606 | Orange <input type="checkbox"/> |
| LN-72042597 | Green <input type="checkbox"/> | LN-72042609 | Assorted* <input type="checkbox"/> |
| LN-72042600 | Pink <input type="checkbox"/> | | |

*Assorted Colors (Blue, Green, Pink, Yellow, Orange)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|------------------|
| 50 | 4 | 30 | | |
| 15 | 12 | 17.5 | | |
| 0.5 | 32 | 8.2 | 6.9 x 3.7 x 2.1 | 17.4 x 9.5 x 5.2 |
| 1.5/2.0 | 32 | 11.6 | | |

Transformer Cube Tube Rack

Low temp storage of large samples,
15 mL & 50 mL tubes

- Quick change over with 2 snap-in inserts that hold 15 and 50 mL tubes
- Rack keyed and numeric wells for tube identification
- Clear Lid included for viewable, stackable storage
- Open space available under insert for extra storage



Pack of 5

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H Dia. | L x W x H cm |
|-------------|--------|-------------|-------|-----------------|---------------------|--------------------|
| LN-72042042 | Purple | 15 | 16 | 18 | 5.25 x 5.25 x 5.125 | 13.2 x 13.2 x 12.9 |
| | | 50 | 9 | 31 | | |

OneRack® Multi

Exclusive design, optimal functionality

The most flexible tube rack designed for use with a
multitude of tube diameters and heights

Holds virtually any brand tube 10-30 mm of
varying heights


Multiple Configurations for Workflow Efficiency

- Interchangeable tube rack inserts hold multiple tube sizes firmly in place when alternating between sample vessels
- Inserts can be pivoted in either an upper or lower position to securely accommodate the size and height of tube being used

Multi-Functional Design

- Open design allows for better filling visibility of tube graduations when filling
- Rack wells have a unique teardrop design that enables better retention of tubes over traditional designs
- Rotate tube rack insert to an angle position to facilitate more comfortable pipetting


Unique Well Design
Tear drop shape allows
for better retention of
tubes



Lower position is ideal height
when working with shorter
microcentrifuge tubes

Alphanumeric grid
for quick sample identification
and retrieval


Ergonomic handles
for ease of transporting
the rack




**Angled position for
more comfortable pipetting**



Included:
1 Base + 6 Tube
Rack Inserts



Two 10-13 mm Tube Racks Fits
5 mL, 7.5 mL, 10 mL Tubes



Two 13-17 mm Tube Racks
Fits 5 mL, 15 mL Tubes



Two 25-30 mm Tube Racks
Fits 25 mL, 50 mL Tubes

OneRack® Multi

US Patent No. 9776188

| Item No. | Color | L x W x H in | L x W x H cm | # Inserts Included | Well Dia./Insert | Wells/Insert |
|-------------|----------------|-----------------|-------------------|-----------------------|------------------|--------------|
| LN-72042891 | Blue/Natural | 12.2 x 4.2 x 3 | 30.9 x 10.7 x 7.7 | 6 (2 ea/ Diameter) | 13 mm | 20 |
| | | | | | 17 mm | 20 |
| | | | | | 30 mm | 6 |
| LN-72042893 | Purple/Natural | 12.2 x 4.2 x 3 | 30.9 x 10.7 x 7.7 | 6 (2 ea/ Diameter) | 13 mm | 20 |
| | | | | | 17 mm | 20 |
| | | | | | 30 mm | 6 |

Adapt-a-Rack® Adaptable Multi-tube Rack

Each well has fingers that flex that adapt to accommodate a variety of tubes sizes while holding them in a firm upright position

- Made of POM
- Racks can be linked together for additional capacity
- Open sides allow for label viewing
- Flat base to hold tubes upright
- Available in solid or a combination of colors for easy coding identification



Pack of 2

| Item No. | Rack #1 Color | Rack #2 Color |
|-------------|---------------|---------------|
| LN-72041316 | Blue/Green | Green/Blue |
| LN-72041319 | Blue/Yellow | Yellow/Blue |
| LN-72041322 | Pink/White | White/Pink |
| LN-72041325 | Blue/Blue | Blue/Blue |
| LN-72041328 | White/White | White/White |
| LN-72041331 | Pink/Pink | Pink/Pink |
| LN-72041334 | Yellow/Yellow | Yellow/Yellow |
| LN-72041337 | Green/Green | Green/Green |

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|------------------|
| 7-50 | 4 | 12-30 | 7.1 x 2.2 x 3 | 18.1 x 5.6 x 7.6 |



Adjustable

Twist the knob to adjust the well sizes for your tubes



Universal Tube Rack

Ideal for holding a wide range of tubes

- Made of polypropylene
- Holds up to 12 similar size tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Side handles for easy transport

Fits
50 mL
Conical
Tubes

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|--------------|-------------|-------|-----------------|-----------------|-----------------|
| LN-72042282 | White/Purple | 15-50 | 12 | 16-30 | 7.6 x 5.5 x 2 | 19.4 x 14 x 5.1 |

Pop-Up Tube Racks

Collapse for compact storage and “pop-up” for quick use

Designed to save space!

- Easy to open, holds 15 and 50 mL tubes
- Easy to close, stackable for compact storage
- Alphanumeric well identification for quick location and retrieval of samples

Pack of 2

| Item No. | Color | Tubes mL | Wells | Well Dia. mm |
|-------------|--------|-------------|--------|-----------------|
| LN-72042264 | Blue | 15, 50 | 21, 12 | 17, 30 |
| LN-72042267 | Pink | | | |
| LN-72042258 | Purple | 15 | 45 | 17 |
| LN-72042261 | Green | 50 | 18 | 30 |

| L x W x H (Expanded) in | L x W x H (Collapsed) in | L x W x H (Expanded) cm | L x W x H (Collapsed) cm |
|----------------------------|-----------------------------|----------------------------|-----------------------------|
| 10 x 5.4 x 2.8 | 10 x 5.4 x 0.8 | 25.5 x 13.7 x 7.2 | 25.5 x 13.7 x 2.1 |



Save Space

Rack easily “Pops” open for easy use



Fold and Snap Tube Rack


Economical durable rack

- Polypropylene racks ship flat in one piece
- Easy-to-assemble, just fold and snap together
- Stable design can be used in water baths
- Open design allows samples to be incubated at the same temperature
- Color code for easy identification



Easy to assemble


Just fold and snap together!



90-Well



60-Well



40-Well

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|--------|-------------|-------|-----------------|-----------------|-------------------|
| LN-72042222 | Blue | | | | | |
| LN-72042228 | White | | | | | |
| LN-72042225 | Red | 5–10 | 90 | 13 | 9.6 x 4.1 x 2.4 | 24.5 x 10.5 x 6.3 |
| LN-72042231 | Yellow | | | | | |
| LN-72042234 | Blue | | | | | |
| LN-72042240 | White | | | | | |
| LN-72042237 | Red | 5–15 | 60 | 16 | 9.6 x 4.2 x 2.8 | 24.5 x 10.8 x 7.2 |
| LN-72042243 | Yellow | | | | | |
| LN-72042210 | Blue | | | | | |
| LN-72042216 | White | | | | | |
| LN-72042213 | Red | 5–30 | 40 | 20 | 9.8 x 4.2 x 2.8 | 24.9 x 10.8 x 7.2 |
| LN-72042219 | Yellow | | | | | |



OneRack® Tube Rack

Tough and highly resistant to acids, bases, solvents and heat

Durable racks are molded in a single, continuous piece so no assembly required and no detachable pieces

Intuitive Design

- Unique tear drop shaped openings offer better retention of tubes
- Molded in coordinates for tube orientation
- Flat labeling area for easier identification
- Stable design for use with water baths

Available in two materials:

Fiberglass-Reinforced Polyoxymethylene (POM)

- Available in half and full rack sizes for 5 different tube diameters
- Offers the highest chemical and heat resistance (-40 to +121°C)

Polypropylene

- Available in full rack sizes for 2 different tube diameters
- Offers good chemical and heat resistance (0 to 110°C)

OneRack® Full-Size POM



| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in (cm) |
|-------------|---------|-------------|-------|-----------------|----------------------|
| LN-72042375 | White | | | | |
| LN-72042378 | Blue | | | | 7.9 x 4 x 2.3 |
| LN-72042381 | Green | 5-10 | 72 | 13 | (20.2 x 10.4 x 5.9) |
| HS27511D | Magenta | | | | |
| LN-72041358 | Yellow | | | | |
| HS27512A | White | | | | |
| LN-72042390 | Blue | | | | 9.8 x 5 x 2.7 |
| LN-72042393 | Green | 5-10 | 72 | 16 | (25 x 12.7 x 7) |
| HS27512D | Magenta | | | | |
| LN-72041361 | Yellow | | | | |
| HS27513A | White | | | | 9.9 x 3.9 x 3.2 |
| LN-72042402 | Blue | | | | (25.2 x 10 x 8.3) |
| LN-72042405 | Green | 10-18 | 40 | 20 | |
| HS27513D | Magenta | | | | |
| HS27514A | White | | | | |
| LN-72042414 | Blue | | | | 11.8 x 4.7 x 3.6 |
| HS27514C | Green | 10-18 | 40 | 25 | (30 x 12 x 9.2) |
| HS27514D | Magenta | | | | |
| HS27515A | White | | | | |
| LN-72042426 | Blue | | | | 11.1 x 4.3 x 3.3 |
| LN-72042429 | Green | 15-50 | 24 | 30 | (28.2 x 11 x 8.5) |
| LN-72042432 | Magenta | | | | |
| LN-72041364 | Yellow | | | | |

OneRack® Half-Size POM



| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in (cm) |
|-------------|---------|-------------|-------|-----------------|----------------------|
| LN-72042327 | White | | | | |
| LN-72042330 | Blue | | | | 4.1 x 4.1 x 2.3 |
| LN-72042321 | Green | 5-10 | 36 | 13 | (10.4 x 10.4 x 5.9) |
| LN-72042336 | Magenta | | | | |
| LN-72042327 | White | | | | |
| LN-72042330 | Blue | | | | 5 x 5 x 2.7 |
| LN-72042333 | Green | 5-10 | 36 | 16 | (12.7 x 12.7 x 7) |
| LN-72042336 | Magenta | | | | |
| LN-72042337 | Yellow | | | | |
| LN-72042339 | White | | | | 4 x 5 x 3.2 |
| HS27503B | Blue | 10-18 | 20 | 20 | (10 x 12.7 x 8.3) |
| LN-72042345 | Green | | | | |
| HS27504A | White | 10-18 | 16 | 25 | 4.7 x 4.8 x 3.6 |
| | | | | | (12 x 12.2 x 9.2) |
| LN-72042363 | White | | | | |
| HS27505B | Blue | | | | 4.3 x 4.3 x 3.3 |
| LN-72042369 | Green | 15-50 | 9 | 30 | (11 x 11 x 8.5) |
| HS27505D | Magenta | | | | |

OneRack® Full-Size Polypropylene


| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in (cm) |
|-------------|---------|-------------|-------|-----------------|----------------------|
| LN-72042435 | White | | | | |
| LN-72042438 | Blue | | | | 7.9 x 4 x 5.9 |
| LN-72042441 | Green | 5-10 | 72 | 13 | (20.2 x 10.4 x 5.9) |
| LN-72042444 | Magenta | | | | |
| LN-72042447 | White | | | | |
| LN-72042450 | Blue | | | | 9.8 x 5 x 2.7 |
| LN-72042453 | Green | 5-10 | 72 | 16 | (25 x 12.7 x 7) |
| LN-72042456 | Magenta | | | | |

Mega Rack®

High-capacity storage for your tubes

Fully assembled, polypropylene mega racks are great for high throughput labs requiring large sample storage


- Sturdy, material for easy and secure transport of blood collection and other tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Place different test groupings into one rack eliminating the need for multiple racks



Mega Racks® 10-13 mm - Holds 5-7 mL sample tubes



| Item No. | Color | Wells | Well Dia. mm | L x W x H in | L x W x H cm | Material |
|-------------|--|-------|--------------|------------------|-------------------|----------|
| LN-72042270 | Blue <input checked="" type="checkbox"/> | 216 | 13 | 13.3 x 9.3 x 3.6 | 33.7 x 23.5 x 9.1 | PP |
| LN-72042273 | White <input type="checkbox"/> | | | | | |
| LN-72041283 | White <input type="checkbox"/> | 432 | 13 | 26.5 x 9.3 x 3.5 | 67.3 x 23.5 x 9.0 | PP |
| LN-72041286 | Blue <input checked="" type="checkbox"/> | | | | | |



Mega Racks® 13-16 mm - Holds 10 mL sample tubes



| Item No. | Color | Wells | Well Dia. mm | L x W x H in | L x W x H cm | Material |
|-------------|--|-------|--------------|-------------------|-------------------|----------|
| LN-72042276 | Blue <input checked="" type="checkbox"/> | 120 | 16 | 13.3 x 10.5 x 3.6 | 33.7 x 26.7 x 9.1 | PP |
| LN-72042279 | White <input type="checkbox"/> | | | | | |
| LN-72041289 | White <input type="checkbox"/> | 240 | 16 | 26.5 x 10.5 x 3.5 | 67.3 x 26.7 x 9.0 | PP |
| LN-72041292 | Blue <input checked="" type="checkbox"/> | | | | | |

Mega Racks® 18-20 mm - Holds 15 mL conical/swab tubes and 20 mm round bottom tubes



| Item No. | Color | Wells | Well Dia. mm | L x W x H in | L x W x H cm | Material |
|-------------|--|-------|--------------|-------------------|-------------------|----------|
| LN-72041156 | White <input type="checkbox"/> | 120 | 20 | 14.9 x 11.9 x 3.6 | 37.7 x 30.1 x 9.1 | PP |
| LN-72041159 | Blue <input checked="" type="checkbox"/> | | | | | |



Peg Tube Racks


Strong, tough, versatile racks for bench use or tube storage

Fiberglass-reinforced polypropylene racks are ideal for bench work or tube storage

- Spacious peg design allows easy viewing of tube contents
- Pegs remain straight even after repeated autoclaving
- Embossed alphanumeric identification for quick location and retrieval of tubes



Unique tear drop well design allows for better retention of tubes



Pack of 2

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|--|----------|-------|--------------|-----------------|-------------------|
| LN-72042246 | White <input type="checkbox"/> | 10 | 50 | 14-17 | 7.4 x 4.1 x 2.7 | 18.8 x 10.4 x 6.9 |
| LN-72042249 | Blue <input checked="" type="checkbox"/> | | | | | |
| LN-72042252 | White <input type="checkbox"/> | 7 | 96 | 10-13 | 3.2 x 5 x 2.7 | 8.2 x 12.6 x 6.9 |
| LN-72042255 | Blue <input checked="" type="checkbox"/> | | | | | |


TracRack® Tube Rack

Take what you need and store the rest

Comprised of a storage box and eight removable tube segments allowing you to easily take what you need - an entire segment or individual sample - and store the rest!

Convenient Rack System

- Eight polypropylene segments hold 1.5–2.0 mL cryo and microcentrifuge tubes (64 tubes)
- Wells feature a center depression stabilizing conical tubes and a cryogenic vial lock for one-handed cap removal
- Printed index letters on each segment make for quick sample ID and retrieval




TracRack® Outer Box

- Made from thin film polypropylene material
- Moisture and mold resistant
- 2 in (5.3 cm) boxes fit in standard freezer racks
- Not recommended for autoclaving

TracRack® Tube Rack

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|--------|----------|-------|--------------|-----------------|-------------------|
| LN-72042009 | Purple | 1.5/2.0 | 64 | 12.3 | 5.1 x 5.1 x 1.4 | 12.9 x 12.9 x 3.5 |



Side Opening
enables faster cooling of samples to ensure integrity

Alphanumeric Identification
for quick location and retrieval of tubes

Sliding Lid
rotates 180° to become a convenient stand the rack sits on

Low-Temp PCR Rack: 96-Well

Slide, set, work & store PCR samples

Multi-purpose rack converts from compact workstation to lidded storage solution accommodating individual and strip 0.2 mL tubes

Versatile & Robust Solution

- Polycarbonate lid serves as cover and elevated rack base - just slide open and pivot lid under rack for compact workstation
- Lid can be used in right or left-handed orientation and is detachable
- Stackable saving valuable bench-top and freezer space
- Multiple colors offered for coding options



US Patent Pending, Community Design No: 004101756-0001-0002


Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|-----------|----------|--------|
| LN-72041707 | Assorted* | 120540 | Green |
| LN-72041708 | Blue | 120541 | Purple |

*Assorted colors (Purple, Blue, Green, Orange, Yellow)



| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|----------|-------|--------------|-----------------|------------------|
| 0.2 | 96 | 5.2 | 4.9 x 3.6 x 1.5 | 12.5 x 9.2 x 3.7 |



Alphanumeric Identification
for quick location and retrieval of tubes

Microtube Storage Boxes

- Hold 1.5 and 2.0 mL microcentrifuge tubes and cryogenic vials




Stylish & Identifiable

- Made of translucent, vibrantly colored polypropylene allows for visibility inside box and color coding identification
- Tube capacity is indicated on the outside of the box within a distinctive wave design
- Quick visual identification with uniform wells and alphanumeric location identification along with a frosted labeling area

Robust

- Easy open and secure lid fitted with either friction fit or snap clasp designs
- Stackable, saving valuable space in the lab



50-Well Microtube Storage Boxes

Pack of 5

| Item No. | Color |
|-------------|--|
| LN-72041154 | Natural <input type="checkbox"/> |
| LN-72041157 | Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> |

*Assorted colors (Blue, Green, Purple, Yellow, Orange)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|------------------|
| 1.5/2.0 | 50 | 12.3 | 5.6 x 3.6 x 2.2 | 14.2 x 9.2 x 5.6 |



81-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|----------------------------------|-------------|--|
| LN-72041160 | Natural <input type="checkbox"/> | LN-72041172 | Yellow <input type="checkbox"/> |
| LN-72041163 | Blue <input type="checkbox"/> | LN-72041175 | Orange <input type="checkbox"/> |
| LN-72041166 | Green <input type="checkbox"/> | LN-72041178 | Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> |
| LN-72041169 | Pink <input type="checkbox"/> | LN-72041181 | Black <input type="checkbox"/> |

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|-----------------|
| 1.5/2.0 | 81 | 12.3 | 5.1 x 5.1 x 1.9 | 13 x 13 x 4.7 |



100-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|--|----------|---------------------------------|
| LN-72041184 | Natural <input type="checkbox"/> | 120202 | Green <input type="checkbox"/> |
| LN-72041187 | Black <input type="checkbox"/> | 120203 | Yellow <input type="checkbox"/> |
| LN-72041190 | Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> | 120204 | Orange <input type="checkbox"/> |
| LN-72041346 | Blue <input type="checkbox"/> | | |

*Assorted colors (Purple, Blue, Green, Yellow, Orange, Pink)

| Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|-------|-----------------|-----------------|-------------------|
| 1.5/2.0 | 100 | 12.3 | 5.6 x 6 x 2.2 | 14.1 x 15.1 x 5.7 |

Cryogenic Vial Boxes & Partitions

For storage of small samples and tubes

- Made from cardboard
- Colors allow for easy identification and sorting purposes
- All box footprints are 13.3 x 13.3 cm, select from two different heights



Cardboard Cryogenic Vial Boxes

Pack of 12

| Item No. | Color | Tubes mL | L x W x H in | L x W x H cm |
|-------------|-----------------------------------|-------------|-----------------|-------------------|
| LN-72042468 | White <input type="checkbox"/> | 1.5/2.0 | 5.2 x 5.2 x 2 | 13.3 x 13.3 x 5 |
| LN-72042471 | White <input type="checkbox"/> | 1.5/2.0 | 5.2 x 5.2 x 3 | 13.3 x 13.3 x 7.5 |
| LN-72042474 | Assorted <input type="checkbox"/> | | | |
| LN-72042477 | Blue <input type="checkbox"/> | 1.5/2.0 | 5.2 x 5.2 x 2 | 13.3 x 13.3 x 5 |
| LN-72042480 | Red <input type="checkbox"/> | | | |



Vent Holes

allow faster cooling and removal of LN₂

Partitions

Partitions must be ordered separately

Pack of 12

| Item No. | Array | Tubes mL | Wells | Well Dia. mm |
|-------------|---------|-------------|-------|-----------------|
| LN-72042486 | 10 x 10 | 1.5/2.0 | 100 | 11.5 |
| LN-72042495 | 9 x 9 | 1.5/2.0 | 81 | 12 |
| LN-72042489 | 8 x 8 | 1.5/2.0 | 64 | 15 |
| LN-72042492 | 7 x 7 | 1.5/2.0 | 49 | 17.5 |

Microtube Storage Boxes

Safe microtube storage

Low profile boxes are available in two sizes to accommodate standard 0.5 mL and 0.2 mL microtubes and vessels of a similar size

- Partitions included in the boxes
- Will fit into most standard freezer racking



Pack of 10

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H cm |
|-------------|--------------------------------|-------------|-------|-----------------|-----------------|-------------------|
| LN-72042895 | White <input type="checkbox"/> | 0.2 | 169 | 9.5 | 5.2 x 5.2 x 1.1 | 13.3 x 13.3 x 2.9 |
| LN-72042897 | White <input type="checkbox"/> | 0.5 | 100 | 12 | 5.2 x 5.2 x 1.6 | 13.3 x 13.3 x 4.1 |



15/50 mL Tube Storage Boxes & Partitions

- Write-on surface for descriptions or self-numbering systems
- Partitions available in two sizes



Cardboard Tube Storage Box & Lid

Pack of 10

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--------------------------------|-----------------|--------------------|
| LN-72042498 | White <input type="checkbox"/> | 5.7 x 5.7 x 4.8 | 14.5 x 14.5 x 12.2 |

Partitions

Partitions must be ordered separately

Pack of 10

| Item No. | Array | Tubes mL | Wells | Well Dia. mm |
|-------------|-------|-------------|-------|-----------------|
| LN-72042504 | 6 x 6 | 15 | 36 | 23 |
| LN-72042501 | 4 x 4 | 50 | 16 | 35 |

Cryogenic Vial Tube Storage Box

- Two-piece set-up boxes with finger holds
- Boxes made from paper coated cardboard
- 5 x 5 format for tubes up to 14 mm in diameter
- Partitions are included in each box



Case of 200

| Item No. | Color | Tubes mL | Wells | Well Dia. mm | L x W x H in | L x W x H mm |
|-------------|-------------------------------|-------------|-------|-----------------|-----------------|-----------------|
| LN-72042499 | Blue <input type="checkbox"/> | 2.0 | 25 | 14 | 3.4 x 3.4 x 1.5 | 87x87x38 |



Sliding Lock


provides secure storage with easily identified lock and unlock symbols.

Slide Rail Design

ensures a secure fit allowing easy placement and removal of slides

Centrally Positioned Numbering

enables simple slide organization with easy viewing of location numbers



100-Place Premium Plus Slide Box

Modern design for maximum functionality

Optimal Organization

- Organize by "project" and "date" molded in on lid
- Inventory card index corresponds to slide numbers and is removable for trouble-free writing on flat surface
- Vibrant translucent colors are easy to identify and enable coding options

Functional Design

- Flexible storage, store vertically or horizontally to efficiently use space
- Molded feet allow secure horizontal stacking of multiple units
- Unique side storage can accommodate a pen, desiccant or other small items needed for working
- Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement

Safe & Secure

- EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- Stylish wave surface texture is easy to hold for safe sample transport

100-Place Premium Plus Slide Box *US Patent No.D904026* Pack of 5

| Item No. | Color | Item No. | Color |
|-------------|--------|-------------|-----------|
| LN-72042811 | Blue | LN-72042814 | Red |
| LN-72042812 | Orange | LN-72042815 | Purple |
| LN-72042813 | Green | LN-72042816 | Assorted* |

*Assorted colors (Purple, Blue, Green, Red, Orange)

| Lining | Lock | L x W x H in | L x W x H cm |
|----------|---------------|-----------------|-------------------|
| EVA Foam | Sliding Latch | 8.2 x 7.7 x 1.3 | 20.7 x 19.5 x 3.2 |

Replacement Inventory Card

Pack of 10

| Item No. | L x H in | L x H cm |
|-------------|-------------|-------------|
| LN-72042822 | 8.5 x 11 | 21.6 x 28 |

Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (76 x 26 mm and 3 x 1 in)

Organize & Identify

- Lining and slots are numbered to correspond to the slide inventory sheet
- Organize by "Project" and "Date" with pre-molded labeling area
- Available in multiple colors for project coding

Durable & Safe

- ABS plastic boxes available with either cork or foam lining
- Designed to securely store slides in grooved slots
- Molded feet secure stacking for easy transport and storage



U.S. Patent No. D904,027

25-Place Slide Box

Easily access slides with friction-fit cover



| Item No. | Color | Lining | Lock | L x W x H in | L x W x H cm |
|-------------|--|--------|--------------|-----------------|------------------|
| LN-72041730 | Blue ■ | Cork | Friction Fit | 5.6 x 3.5 x 1.4 | 14.1 x 8.8 x 3.5 |
| LN-72041733 | Green ■ | | | | |
| LN-72041736 | Red ■ | | | | |
| LN-72041739 | White □ | Foam | Friction Fit | 5.6 x 3.5 x 1.4 | 14.1 x 8.8 x 3.5 |
| LN-72041742 | Blue ■ | | | | |
| LN-72041745 | Red ■ | | | | |
| LN-72041748 | White □ | | | | |

50-Place Slide Box

Thumb-latch lock mechanism



| Item No. | Color | Lining | Lock | L x W x H in | L x W x H cm |
|-------------|--|--------|-------------|-----------------|------------------|
| LN-72041799 | Blue ■ | Cork | Thumb-Latch | 8.2 x 3.4 x 1.4 | 20.8 x 8.6 x 3.5 |
| LN-72041802 | Red ■ | | | | |
| LN-72041805 | White □ | | | | |

100-Place Slide Box

Closes securely with a rust-resistant nickel-plated clasp and pin hinges



| Item No. | Color | Lining | Lock | L x W x H in | L x W x H cm |
|-------------|--|--------|---------------------|-----------------|-------------------|
| LN-72041769 | Blue ■ | Cork | Nickel-plated Clasp | 8.2 x 6.9 x 1.3 | 20.8 x 17.5 x 3.4 |
| LN-72041772 | Green ■ | | | | |
| LN-72041775 | Red ■ | | | | |
| LN-72041778 | Yellow ■ | | | | |
| LN-72041781 | White □ | | | | |
| LN-72041784 | Grey ■ | | | | |
| LN-72041787 | Black ■ | Foam | Nickel-plated Clasp | 8.2 x 6.9 x 1.3 | 20.8 x 17.5 x 3.4 |
| LN-72041790 | Blue ■ | | | | |
| LN-72041793 | Red ■ | | | | |
| LN-72041796 | White □ | | | | |

100-Place Tall Slide Box

Additional height for 3 x 2 in slides
Stackable with 100-Place Premium Slide Boxes



| Item No. | Color | Lining | Lock | L x W x H in | L x W x H cm |
|-------------|--|--------|---------------------|-----------------|-----------------|
| LN-72041610 | Blue ■ | Cork | Nickel-plated Clasp | 8.2 x 6.9 x 2.4 | 20.8 x 17.5 x 6 |



Storage Rack for Microscope Slide Boxes

Create a file system for 50-Place and 100-Place slide boxes

Storage rack makes organization of samples easy

- Clear acrylic construction allows the labels to be seen on the slide boxes
- Place boxes in rows or columns by color to create an easy coding system

| Item No. | Color | Capacity | L x W x H in | L x W x H cm |
|-------------|-------|----------|-----------------|--------------------|
| LN-72041727 | Clear | 15 Bins | 28 x 6.7 x 9 | 71.1 x 17.1 x 22.8 |

Standard Microscope Slide Boxes

Constructed of durable polypropylene



Slide boxes hold either 25 or 100 standard slides (76 x 26 mm and 3 x 1 in)

Improved Identification

- Molded slot numbers correspond directly to the slide inventory card
- Organize with "Project" and "Date" pre-molded labeling area

Secure & Stackable

- Hinged lids are held securely closed with a latch lock
- Molded feet secure stacking for easy transport and storage

25-Place Slide Box

| Item No. | Color | Lining | Lock | L x W x H in (cm) |
|-------------|---|--------|------------|-------------------------------------|
| LN-72041751 | Blue | None | Latch Lock | 5.6 x 3.6 x 1.2 (14.1 x 9.2 x 3) |
| LN-72041754 | Red | | | |
| LN-72041757 | White | | | |
| LN-72045307 | 25-Place Replacement Inventory Card, Pack of 10 | | | 5.3 x 3.1 (13.4 x 8) |

100-Place Slide Box

| Item No. | Color | Lining | Lock | L x W x H in (cm) |
|-------------|--|--------|------------|--------------------------------------|
| LN-72041760 | Blue | None | Latch Lock | 8.3 x 6.7 x 1.5 (21 x 16.9 x 3.7) |
| LN-72041763 | Red | | | |
| LN-72041766 | White | | | |
| LN-72045308 | 100-Place Replacement Inventory Card, Pack of 10 | | | 3.9 x 6.1 (10 x 15.5) |

Low Temperature Microscope Slide Boxes

Store in ultra-low temperatures with confidence



Holds standard slides (76 x 26 mm and 3 x 1 in)

Durable Construction

- Made of polycarbonate to store in ultra-low temperatures
- Uniform grooved slide slots separate slides preventing slide contact
- Molded feet enable secure stacking for easy transport and storage

Durable Construction

- Low temperature resistant numbered foam lining corresponds with identification sheet
- Organize with "Project" and "Date" pre-molded labeling area



Low Temperature Microscope Slide Boxes

| Item No. | Color | Capacity | Lock | L x W x H in (cm) |
|-------------|--------|-----------|-----------------------|--|
| LN-72042843 | Blue | 25-Place | Friction Fit | 5.6 x 3.5 x 1.4 (14.1 x 8.8 x 3.5) |
| LN-72041724 | Purple | 100-Place | Stainless Steel Clasp | 8.2 x 6.9 x 1.3 (20.8 x 17.5 x 3.4) |

Cardboard Slide Tray

Cardboard slide trays hold standard slides (76 x 26 mm and 3 x 1 in)

Cardboard slide folders offer a convenient solution for safely working, transporting and short-term storing of slides

Organize & Protect

- Index on the front cover corresponds with the number grid inside for identification
- Recessed slide compartments and inside covers for specimen protection
- Slides cannot be accessed when cover is closed reducing contamination

Sturdy Construction

- Compartment holder is made from a double layer of densely pressed cardboard
- Cover is single layer cardboard
- Front protection flaps are attached with durable tape



| Item No. | Color | Capacity | L x W x H in | L x W x H cm |
|-------------|----------|----------|------------------|-------------------|
| LN-72041726 | Black ■ | | | |
| LN-72042795 | Blue ■ | | | |
| LN-72042798 | Green ■ | 20-Place | 8.2 x 13.4 x 0.4 | 20.8 x 34.1 x 0.9 |
| LN-72041340 | Red ■ | | | |
| LN-72041343 | Yellow ■ | | | |
| LN-72042801 | Black ■ | 30-Place | 8.2 x 12.9 x 0.4 | 20.8 x 32.8 x 1 |

Slide Tray Clasps

Effortlessly locate, retrieve and keep slide folders organized

- Assorted color clasps are ideal for color coding identification and to keep slide folder flaps securely closed
- Clasps can be labeled with a permanent marker or a label can be attached



Slide Tray Clasps

Pack of 100

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--|-----------------|-----------------|
| LN-72042804 | Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> | 2.2 x 1.7 x 0.6 | 5.6 x 4.4 x 1.6 |
| LN-72042804 | Clear <input type="checkbox"/> | | |

*Assorted colors (Clear, Grey, Blue, Green, Red)

Slide Mailers

Safe and secure transport of standard microscope slides (76 x 26 mm and 3 x 1 in)

- Heavy cardboard mailers for safe transport
- Cover that opens for viewing and thumb grooves for easy slide removal
- Securely holds slides with or without cover glass
- Flat surface area is ideal for labeling or writing necessary information



Cardboard Slide Mailers

| Item No. | Color | Capacity | UOM | L x W x H in | L x W x H mm |
|-------------|-----------|----------|------------|-----------------|-----------------|
| LN-72042786 | Natural ■ | 1-Place | Pack of 25 | 3.7 x 1.7 x 0.2 | 95 x 42 x 5 |
| LN-72042789 | Natural ■ | 2-Place | Pack of 36 | 4.1 x 3.1 x 0.2 | 103 x 79 x 5 |
| LN-72042784 | White □ | 4-Place | Pack of 20 | 4.7 x 5.7 x 0.3 | 120 x 140 x 6 |



Cool Containers

Durable containers for chilling temperature-sensitive samples

Suitable for low temperature use. Ideal for ice, ice slurries and CO₂ pellets

- Made from polyurethane with metal inserts for stability
- Close-fit lid keeps evaporation to a minimum and features a comfortable recessed knob
- Tapered design for easy stacking
- Precise pour spout and large handles



Lid nests under bucket
for convenience and to save space



| Item No. | Color | Style | Capacity liters | L x W x H in | L x W x H cm |
|-------------|---------|--------|--------------------|-------------------|---------------------|
| LN-72042579 | Black ■ | Bucket | 4 | 7.1 (H) x 9.2 (D) | 18.1 (H) x 23.4 (D) |
| LN-72042582 | Green ■ | Pan | 1 | 10.5 x 7.7 x 4.5 | 26.7 x 19.6 x 11.4 |
| LN-72042585 | Blue ■ | Pan | 4 | 15 x 10 x 5.5 | 38.1 x 25.4 x 14 |
| LN-72042588 | Red ■ | Pan | 9 | 20.5 x 11.5 x 6.4 | 52.1 x 29.2 x 16.3 |

LevGo® Eco smartFunnel™ Paper Funnel



A clean transfer every time

- Use with virtually any liquid or powder (not recommended for hot liquids)
- Often used to collect sputum samples
- Environmentally friendly paper product that will naturally biodegrade in the standard municipal waste stream
- Made of FSC® Certified, FDA food grade, clay coated heavy grade paper
- Innovative "backrest" allows hands-free use when pouring viscous liquids
- Maximum flexibility - cut stem to any desired opening size or shape
- Apply a minimum of pressure and the funnel will stay in position



Pack of 25

| Item No. | Size | Capacity mL | H x Inlet x Stem in/mm |
|-------------|------|-------------|--|
| LN-72041367 | Mini | 155 | 5 x 21/4 x 2 x 1/4 in ² / 127 x 57 mm ² x 64 mm ² |

Lead Rings

Stabilize vessels and labware on the bench or in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit



| Item No. | Shape | Color | Flask Sizes mL | ID in | ID cm | OD in | OD cm | Weight g |
|-------------|-------|-------|----------------|-------|-------|-------|-------|----------|
| LN-72042759 | C | Blue | 125–500 | 1.9 | 4.2 | 3 | 7.5 | 215 |
| LN-72042771 | O | | | 1.9 | 4.8 | 3 | 7.7 | 210 |
| LN-72042762 | C | Blue | 250–1000 | 2 | 5.1 | 3.5 | 9 | 465 |
| LN-72042774 | O | | | 2 | 5.1 | 3.5 | 9 | 490 |
| LN-72042765 | C | Blue | 500–2000 | 2.1 | 5.4 | 3.9 | 9.8 | 595 |
| LN-72042777 | O | | | 2.2 | 5.7 | 3.9 | 9.8 | 610 |
| LN-72042768 | C | Blue | 1000–4000 | 2.6 | 6.6 | 4.4 | 11.3 | 890 |
| LN-72042780 | O | | | 2.8 | 7 | 4.5 | 11.5 | 935 |

*Weights may vary by +/- 3%

Safety Bottle Carriers

Secure and convenient bottle transport

- Ergonomic handle prevents hand fatigue and related accidents when carrying heavier loads
- Handle swings down to either side for easier loading and unloading
- Centering ring stabilizes bottle in the center of carrier for safer, steadier transport
- Carriers come in four colors with designated labeling areas for easy identification
- Made of chemical resistant thermoplastic elastomer - PP blend



Available in 4 Colors
for easy identification and organization

Safety Bottle Carrier, 4.5 L

| Item No. | Color | Dimensions (Dia. x H) in | Dimensions (Dia. x H) cm |
|-------------|-------|---|--|
| LN-72041378 | Red | 7.2 x 10.5 (16 in height, handle extended) | 18.3 x 26.7 (40.6 cm height, handle extended) |
| LN-72041379 | Black | | |
| LN-72041380 | Blue | | |
| LN-72041381 | Green | | |

Safety Bottle Carrier, 1.5 L

| Item No. | Color | Dimensions (Dia. x H) in | Dimensions (Dia. x H) cm |
|-------------|-------|--|--|
| LN-72041382 | Red | 4.5 x 6.5 (10 in height, handle extended) | 11.4 x 16.5 (25.4 cm height, handle extended) |
| LN-72041383 | Black | | |
| LN-72041384 | Blue | | |
| LN-72041385 | Green | | |



Movable Handle
swings down to allow
clear access to vessel



Biohazard Disposal Bags

For personal or general lab use, sizes to suit all needs

- Polypropylene
- A range of biohazard collection bags printed with biohazard warning symbol and precautions in English, Spanish and French
- Includes an indicator patch that darkens when autoclaved



Biohazard Disposal Bags

| Item No. | Color | Thickness mil | Capacity | Packaging | L x W in | L x W cm |
|-------------|-------|------------------|----------|-----------|-------------|-------------|
| LN-72041055 | Red ■ | 1.6 | 2 L | 500 | 8 x 12 | 20.3 x 30.5 |
| LN-72041058 | | | 9 L | 200 | 14 x 19 | 35.6 x 48.3 |
| LN-72041061 | | | 19 L | 200 | 19 x 23 | 48.3 x 58.4 |
| LN-72041064 | | | 40 L | 200 | 25 x 33 | 63.5 x 83.8 |

HDPE Coated Steel Wire Bag Holder

| Item No. | L x W x H in | L x W x H cm |
|-------------|-----------------|-------------------|
| LN-72041034 | 3.6 x 5 x 8.4 | 9.1 x 12.8 x 21.3 |

*Holds Biohazard Disposal Bag HS10320



Epoxy Coated
HDPE Wire Holder



Autoclavable Disposal Bags

Clear bags imprinted with blue biohazard labeling designed for safe disposal

- Constructed of 2 mil thick polypropylene, allowing for repeated autoclaving
- Front write-on panel
- 4 inner packs of 50 bags
- Strong side seams for securely holding items



Case of 200

| Item No. | Color | Capacity | L x W in | L x W cm |
|-------------|---------|----------|-------------|-------------|
| LN-72041022 | Clear □ | 14 L | 12.2 x 26 | 31 x 66 |
| LN-72041025 | | 26 L | 16.3 x 23.6 | 41.5 x 60 |
| LN-72041028 | | 69 L | 24 x 31.9 | 62 x 81 |

Sharps Containers

Easy to use and transport

- Clearly labeled with the Biohazard warning on one side
- Provides safe disposal of syringes, needles and other small sharps
- Stackable to reduce storage space before use
- 1 L Capacity (quart) with a clear lid to view fill level
- Slide access lid provides temporary and final closures for safety and security
- Built-in needle removal port, prevents direct contact with sharps
- Compatible with the Droplet blood collection tray, Totemax® and blood collection tray

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|----------|-----------------|--------------------|
| LN-72041307 | Red ■ | 4.5 x 4.5 x 7.5 | 11.5 x 11.5 x 19.2 |
| LN-72041310 | Yellow ■ | | |



PPE Dispensers & Holders



Soft Cover Dispenser

Store and dispense with ease

Translucent acrylic dispensers provide easy access to safety accessories

- Hinged lid protects contents from dust and debris
- Set on counter or wall mount (screws included)
- Great for holding other personal items



| Item No. | Color | L x W x D in | L x W x D cm |
|-------------|--------------------------------|-----------------|--------------------|
| LN-72041088 | Clear <input type="checkbox"/> | 8 x 8 x 16 | 20.3 x 20.3 x 40.6 |

Earplug Dispenser

Translucent acrylic dispensers provide easy access to safety accessories

- Holds approximately 200 pairs of foam earplugs
- Set on counter or wall mount (screws included)



| Item No. | Color | L x W x D in | L x W x D cm |
|-------------|-------------------------------------|-----------------|--------------------|
| LN-72041085 | Neon Green <input type="checkbox"/> | 8 x 8 x 16 | 20.3 x 20.3 x 40.6 |

3-in-1 Dispenser

Three separate lift-up lids make loading and dispensing supplies easy

Three-compartment dispenser can be used to organize safety glasses, earplugs, soft covers (including bouffant caps, beard covers, shoe covers) or any small items that need to be kept within reach

The unit includes screws for wall mounting or can be set on a counter

- Middle compartment has reach-in hole for dispensing
- Durable acrylic is see-through so product levels are always visible

3
Compartment



| Item No. | Color | L x W x D in | L x W x D cm |
|-------------|-------------------------------|-----------------|--------------------|
| LN-72041094 | Grey <input type="checkbox"/> | 16 x 16 x 8 | 40.6 x 40.6 x 20.3 |

Safety Glasses Dispensers

Ergonomically angled for easy dispensing

Translucent acrylic dispensers provide easy access to safety accessories

- Hinged lid protects glasses from dust and debris
- Set on counter or wall mount (screws included)
- Large dispenser holds approx. 20 glasses, small approx. 10 glasses
- Great for holding other personal items

| Item No. | Color | L x W x D in | L x W x D cm |
|-------------|------------------------------------|-----------------|--------------------|
| LN-72041091 | Neon Blue <input type="checkbox"/> | 8 x 6 x 6 | 20.3 x 15.2 x 15.2 |
| LN-72041082 | Neon Blue <input type="checkbox"/> | 8 x 8 x 16 | 20.3 x 20.3 x 40.6 |



PPE Dispensers & Holders

Workstation Storage Bin

The essential combination accessory holder

- Organize your most commonly used laboratory disposables in one central location
- Acrylic unit has room to store gloves, Kimwipes® and Parafilm® M
- Use bins to store tubes, tips, etc.

| Item No. | Color | L x W x H in | L x W x H cm |
|-------------|--------------------------------|-------------------|------------------|
| LN-72042093 | Clear <input type="checkbox"/> | 10.4 x 5.4 x 17.3 | 26.5 x 13.7 x 44 |

*Parafilm® is a registered trademark of Amcor Limited
Kimwipes® is a registered trademark of Kimberly-Clark*



Lab Wipe Holders

Holds Kimwipes® and other major brands

Clear acrylic holders fit boxes of small or large lab wipes

- The smaller version attaches with double-sided adhesive tape and the larger version with screws
- Mount on a shelf or wall
- Mounting materials included

| Item No. | Description | Color | L x W x H in (cm) |
|-------------|--------------|--------------------------------|--|
| LN-72042087 | Small Holder | Clear <input type="checkbox"/> | 4.6 x 4.8 x 3.3 (11.8 x 12.2 x 8.4) |
| LN-72042090 | Large Holder | Clear <input type="checkbox"/> | 12.6 x 5 x 3.8 (32 x 12.7 x 9.7) |

Kimwipes® is a registered trademark of Kimberly-Clark

Benchtop Dispensing Bins

Work great with gloved hands, bright colors

Ergonomically designed with large, angled openings and oversized lip

- Lid remains open when tilted back
- Nonskid rubber feet and pre-drilled holes for wall mounting
- Green bin includes a removable divider
- Mounting screws included

| Item No. | Color | Compartments | L x W x H in | L x W x H cm |
|-------------|---|--------------|-----------------|------------------|
| LN-72042072 | Neon Green ■ | 2 | 7.1 x 6.1 x 6.5 | 18 x 15.5 x 16.5 |
| LN-72042075 | Neon Blue ■ | 1 | 6.1 x 6.1 x 6.7 | 15.5 x 15.5 x 17 |



Clear Acrylic Glove Box Holders

Keep gloves organized & ready for immediate use

- Free up bench-top space by mounting holders on wall for easy dispensing of gloves
- Universal fit for most major brands of glove boxes
- Mounting screws included
- Includes a polyester foam insert to securely hold a range of glove box sizes—snug fit ensures easy removal of gloves from the box

| Item No. | Capacity | L x W x H in | L x W x H cm |
|-------------|----------|-------------------|--------------------|
| LN-72042171 | 1 | 10.3 x 5.9 x 4.5 | 26.3 x 14.6 x 11.4 |
| LN-72042174 | 2 | 11.4 x 10.3 x 4.5 | 29.2 x 26 x 11.4 |
| LN-72042177 | 3 | 11.6 x 15.2 x 4.3 | 29.4 x 38.8 x 11.6 |



HDPE Coated Wire Glove Box Holders

Keep gloves organized & ready for immediate use

- HDPE-coated steel wire racks are an economical solution
- Free up bench-top space by mounting holders on wall for easy dispensing of gloves
- Universal fit for most major brands of glove boxes
- Mounting screws included

| Item No. | Capacity | L x W x H in | L x W x H cm |
|-------------|----------|------------------|------------------|
| LN-72042126 | 1 | 5.5 x 8.2 x 4.3 | 14 x 20.9 x 10.8 |
| LN-72042129 | 3 | 9.8 x 17.9 x 4.3 | 25 x 45.4 x 10.8 |

Modular Glove Box Holders

Available infused with Microban® for anti-microbial protection

The anti-microbial version (blue) is ideal for use in clean rooms, hospitals, and microbiological work areas. The additive prevents growth of molds and bacteria

- High-impact, durable resin is resistant to harsh cleaning solutions
- Universal fit for most major brands of glove boxes
- Modular design allows single boxes or horizontal stacking for multiple box holders
- Mounting screws included

| Item No. | Description | Color | US Patent No. 6669157 Pack of 3 | |
|-------------|----------------|--|---------------------------------|-------------------|
| | | | L x W x H in | L x W x H cm |
| LN-72041445 | Standard | White <input type="checkbox"/> | 10 x 6.1 x 3.8 | 25.4 x 15.6 x 9.7 |
| LN-72042180 | Anti-microbial | Blue <input checked="" type="checkbox"/> | | |



Heat Sealer with Cutter

Quickly and efficiently seal polyethylene or polypropylene bags, sheets or flat tubing from 0.01 to 0.8 mm thick


Sliding blade removes excess material 19 mm from the seal

- Can seal a bag up to 12 in (30.5 cm) long
- Heating element measures 31 cm long and is only hot when sealing
- Includes one replacement part kit
- Additional replacement part kits can be ordered separately
- For use with PE and PP material items only. Not compatible with combination plastic/paper pouches



Replacement Part Kit:

- 2 extra blades,
- 2 extra heating elements,
- 2 extra PTFE cloths



| Item No. | Power Supply | L x W X H in | L x W X H cm |
|-------------|----------------------------------|------------------|-------------------|
| LN-72045632 | 220-230 VAC, Cont. Europe/Israel | 17.6 x 3.4 x 6.5 | 44.6 x 8.6 x 16.6 |
| LN-72045640 | Replacement Part Kit | | |

Parafilm® M Sealing Film

Protect the integrity of your laboratory and research applications

Parafilm® M is a semi-transparent, flexible, thermoplastic film with two-way stretch

- Protect the integrity of research applications with a tight moisture-barrier that holds a strong seal and doesn't release from glassware, despite off-gassing
- Reduce the risk of tampered data with a covering that is odorless, colorless and semi-transparent
- Prevents moisture and volume loss due to material composition and tight seal
- Flexible and self-sealing material allows it to cling around irregular shapes and surfaces
- Stretches more than 200% of original length
- Resistant up to 48 hours against many polar substances such as saline solutions, inorganic acids and alkaline solutions
- Two color options available in 2" rolls (natural and purple) for color coding purposes



| Item No. | Color | | Film Width | | Roll Length | |
|-------------|---------|-------------------------------------|------------|------|-------------|------|
| | | | in | cm | ft | m |
| LN-72042111 | Natural | <input type="checkbox"/> | 2 | 5 | 250 | 76.2 |
| LN-72045320 | Purple | <input checked="" type="checkbox"/> | 2 | 5 | 250 | 76.2 |
| BR-70701605 | Natural | <input type="checkbox"/> | 4 | 10 | 125 | 38.1 |
| BR-70701606 | Natural | <input type="checkbox"/> | 4 | 10 | 250 | 76.2 |
| BR-70701501 | Natural | <input type="checkbox"/> | 20 | 50.8 | 50 | 15.2 |

Parafilm® M is a registered trademark of Amcor Limited

Box Top Dispenser for Parafilm® M Sealing Film

Easy application

Simply place the dispenser on top of the Parafilm® M box and pull the sealing film through, there is no need to remove rolls from the packaging

- Accommodates two 50 mm or one 100 mm Parafilm® M rolls
- Built in safety blade slides and cuts with no danger to fingertips and eliminates the need for scissors
- Supplied with two safety blades, additional replacement blades are available for purchase
- Made with ABS plastic making the unit easy to clean

| Item No. | | L x W x H in | L x W x H cm |
|-------------|-----------|-----------------|-------------------|
| LN-72045155 | Blue ■ | | |
| LN-72045156 | Purple ■ | 4.9 x 6.8 x 3.9 | 12.3 x 17.2 x 9.9 |
| LN-72045318 | Natural □ | | |

| Item No. | |
|-------------|---|
| LN-72041582 | Replacement Blades for Box Top Dispenser, Pack of 3 |

Parafilm® M is a registered trademark of Amcor Limited



ABS Dispenser for Parafilm® M Sealing Film

Store, dispense, and cut Parafilm® M, tape, labels, and Tough-Tags®

- ABS plastic dispenser designed with integrated safety razor that slides smoothly creating a clean cut
- Supplied with two safety blades, additional replacement blades are available for purchase
- Dispenser holds one or two 50 mm (2-inch) rolls or one 100 mm (4-inch) roll

| Item No. | | L x W x H in | L x W x H cm |
|-------------|---------|-----------------|------------------|
| LN-72042099 | Blue ■ | | |
| LN-72042102 | Green ■ | 4.7 x 6.1 x 6.7 | 12 x 15.6 x 17.1 |

| Item No. | |
|-------------|---|
| LN-72042108 | Replacement Blades for ABS Dispenser, Pack of 5 |

Parafilm® M is a registered trademark of Amcor Limited
Tough-Tags® is a registered trademark of Diversified Biotech

Acrylic Dispenser for Parafilm® M Sealing Film

Store, dispense, and cut Parafilm® M, tape, labels, and Tough-Tags®

- Strong, crystal-clear acrylic dispenser makes your film supply visible and features a serrated cutting blade
- Rubber feet secures in place while cutting
- Dispenser holds one or two 50 mm (2-inch) rolls or one 100 mm (4-inch) roll

| Item No. | L x W x H in | L x W x H cm |
|-------------|-----------------|--------------------|
| LN-72042096 | 4.7 x 6.6 x 6.9 | 11.9 x 16.8 x 17.5 |

Parafilm® M is a registered trademark of Amcor Limited
Tough-Tags® is a registered trademark of Diversified Biotech



Partner per l'Italia



STEDA SRL
Marketing Trading Logistics

Steda srl
Via Cavour 7 , 20814 VAREDO (MB)
www.stedalab.it