




Slide Rail Design
ensures a secure fit allowing easy placement and removal of slides

Centrally Positioned Numbering
enables simple slide organization with easy viewing of location numbers

Sliding Lock
provides secure storage with easily identified lock and unlock symbols.



100-Place Premium Plus Slide Box

Modern design for maximum functionality

Optimal Organization

- Organize by "project" and "date" molded in on lid
- Inventory card index corresponds to slide numbers and is removable for trouble-free writing on flat surface
- Vibrant translucent colors are easy to identify and enable coding options

Functional Design

- Flexible storage, store vertically or horizontally to efficiently use space
- Molded feet allow secure horizontal stacking of multiple units
- Unique side storage can accommodate a pen, desiccant or other small items needed for working
- Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement

Safe & Secure

- EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- Stylish wave surface texture is easy to hold for safe sample transport

100-Place Premium Plus Slide Box

Pack of 5

Item No.	Color	Item No.	Color
120576	Blue	120579	Red
120577	Orange	120592	Purple
120578	Green	120593	Assorted*

*Assorted colors (Purple, Blue, Green, Red, Orange)

Lining	Lock	L x W x H in	L x W x H cm
EVA Foam	Sliding Latch	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2

Replacement Inventory Card

Pack of 10

Item No.	L x H in	L x H cm
120618	8.5 x 11	21.6 x 28

Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (76 x 26 mm and 3 x 1 in)

Organize & Identify

- Lining and slots are numbered to correspond to the slide inventory sheet
- Organize by "Project" and "Date" with pre-molded labeling area
- Available in multiple colors for project coding

Durable & Safe

- ABS plastic boxes available with either cork or foam lining
- Designed to securely store slides in grooved slots
- Molded feet secure stacking for easy transport and storage



25-Place Slide Box

Easily access slides with friction-fit cover



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15989A	Blue <input checked="" type="checkbox"/>	Cork	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989B	Green <input checked="" type="checkbox"/>				
HS15989C	Red <input checked="" type="checkbox"/>				
HS15989E	White <input type="checkbox"/>	Foam	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989M	Blue <input checked="" type="checkbox"/>				
HS15989P	Red <input checked="" type="checkbox"/>				
HS15989R	White <input type="checkbox"/>				

50-Place Slide Box

Thumb-latch lock mechanism



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15996A	Blue <input checked="" type="checkbox"/>	Cork	Thumb-Latch	8.2 x 3.4 x 1.4	20.8 x 8.6 x 3.5
HS15996B	Red <input checked="" type="checkbox"/>				
HS15996C	White <input type="checkbox"/>				

100-Place Slide Box

Closes securely with a rust-resistant nickel-plated clasp and pin hinges



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15994A	Blue <input checked="" type="checkbox"/>	Cork	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994B	Green <input checked="" type="checkbox"/>				
HS15994C	Red <input checked="" type="checkbox"/>				
HS15994D	Yellow <input checked="" type="checkbox"/>				
HS15994E	White <input type="checkbox"/>				
HS15994F	Grey <input type="checkbox"/>				
HS15994G	Black <input checked="" type="checkbox"/>	Foam	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994M	Blue <input checked="" type="checkbox"/>				
HS15994P	Red <input checked="" type="checkbox"/>				
HS15994R	White <input type="checkbox"/>				

100-Place Tall Slide Box

Additional height for 3 x 2 in slides
Stackable with 100-Place Premium Slide Boxes



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15979A	Blue <input checked="" type="checkbox"/>	Cork	Nickel-plated Clasp	8.2 x 6.9 x 2.4	20.8 x 17.5 x 6